

Shiley™

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Shiley™

Oral/Nasal Tracheal Tube Cuffless Non DEHP, Murphy Eye

Instructions for Use

Description

The Shiley™ oral/nasal cuffless non DEHP pediatric endotracheal tube with Murphy eye is a translucent tube with a hooded Murphy eye tip. The tube incorporates a Magill curve and features a radiopaque line. The tube is supplied sterile with a standard 15 mm connector. In addition to the standard depth markings (length markings) indicating the distance to the distal tip in centimeters, a second series of reference lines, glottic depth marks, are printed at the distal end. The glottic depth marks aid tube positioning during intubation and their distance from the distal tip varies with tube size as summarized in the table on page 2.

This product and its packaging are not made with natural rubber latex or DEHP.

Sterilization is by ethylene oxide.

Indications/Intended Use/Purpose

The oral/nasal cuffless non DEHP pediatric endotracheal tube with Murphy eye is intended for oral or nasal intubation of the trachea and is indicated for airway management.

Contraindications

None known.

WARNINGS

- This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.
- If the device is damaged or the integrity of the sterile barrier has been compromised, do not use the product and contact your Covidien representative for further information.
- This device must only be used by clinicians trained in its use.
- Avoid contact of the tube with laser beams during laser surgery or electrosurgical active electrodes. Contact can result in tube ignition, combustion, airway fire and patient injury. Consider using a laser-resistant tracheal tube if airway laser surgery is planned.
- If the tracheal tube is lubricated prior to insertion, follow manufacturer's lubricant application instructions. Avoid lubricant entering and occluding the tube lumen.
- It is essential to verify that the tube position remains correct after intubation, especially when a patient's position or the tube placement is altered. Correct any tube mal-position immediately.
- When properly sized, there may be a slight leak around the tracheal tube resulting in a loss of tidal volume with positive pressure ventilation. Compensatory adjustments may need to be made as dictated by expert clinical judgment.
- Use of a reinforced tracheal tube should be considered if extreme flexion of the neck (chin-to-chest), movement of the patient (e.g., to a lateral or prone position), or tube compression is anticipated after intubation.

CAUTIONS

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- The tube is for single patient use. Duration should not exceed twenty-nine (29) days. Covidien has not substantiated use of these devices beyond 29 days. Decisions about tracheal tube changes should be made by the responsible physician or designate using accepted medical techniques and judgement.
- The user should be alert for anatomical variations, and reliance on the reference lines should not be substituted for expert clinical judgment.
- The 15 mm connector is seated so that it can be removed with effort if cutting of the tube is desired. Follow the Instructions if cutting is considered. Do not use lubricating jelly to ease connector reinsertion, as it may contribute to accidental disconnection.
- Use only with equipment having standard 15 mm fittings.
- During storage, avoid sunlight.
- The device is MR safe.

Instructions

1. Use medical and clinical judgment for selection of the appropriate type and size endotracheal tube for each individual patient.
2. Visually inspect package integrity and device prior to use with consideration given to the specific WARNINGS and CAUTIONS stated in this insert. Remove the endotracheal tube from its protective package.
3. If shortening of the endotracheal tube is needed, use medical judgment to determine the appropriate tube length for each individual patient. Before intubating the patient, cut the tube, reinsert the 15 mm connector, and verify the connector is firmly seated in the tracheal tube.
4. If using a stylet to facilitate intubation, verify the stylet does not extend beyond the end of the tube and can be easily removed before intubation. Take care not to damage the stylet sheath during reshaping, insertion or removal. If the stylet sheath is damaged, do not use.
5. Intubate the patient using currently accepted medical techniques with consideration given to the specific WARNINGS and CAUTIONS stated in this product insert.
6. Attach the 15 mm connector to the breathing circuit. Ensure the connector is firmly seated in the breathing circuit to prevent disconnection during use. Consider use of a swivel adapter with ventilator breathing circuits.
7. Follow currently accepted medical techniques to confirm that the endotracheal tube has been placed into the trachea, and not inadvertently placed into the esophagus or a mainstem bronchus.
8. After intubation and confirmation of correct placement, secure the tracheal tube following currently accepted medical techniques.
9. To remove secretions, follow currently accepted medical techniques to suction the endotracheal tube lumen.
10. Extubate the patient following currently accepted medical techniques.
11. Dispose of the tracheal tube in accordance with current medical standards and applicable national regulations for biologically hazardous waste.

Adverse Events

Adverse events associated with the procedural placement and use of this endotracheal tube in descending order of severity and frequency are: respiratory failure or respiratory distress; tissue damage or trauma; dental injury; post-operative sore throat; pain/discomfort; reintubation; pneumothorax; anaphylactic reaction; aspiration; infection; delay of treatment. Any serious incident that occurs in relation to the device or its accessories should be reported to Covidien and to the national competent authority.

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	Sterilized using ethylene oxide
	Single use
	Do not use if package is damaged
	Consult instructions for use
	Caution
	Keep away from sunlight
	Authorized representative in the European Community
	Manufacturer
	Use-by date
	Lot number
	Date of manufacture

Page Left Intentionally Blank

Page Left Intentionally Blank

Rx
ONLY

CE
0123

Part No. PT00100292 Rev A 2019-11

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG.

TM* brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

© 2019 Covidien.

 Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.

 Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands.

www.covidien.com