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

en 3M™ Ranger™ Fluid Warming Irrigation Set

Indications for Use

The 3M™ Ranger™ Irrigation Warming System is intended to warm irrigation fluids.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas where irrigation fluid warming is required.

Explanation of Signal Word Consequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

WARNING: To reduce the risks of biohazard exposure, cross-contamination or infection:

- Do not reuse. This product is designed for single patient use only.
- Do not use for IV infusion; use for irrigation only
- Do not exceed 245 mm/Hg pressure; ensure inline hand pressure pumps or other pumps do not exceed 245 mm/Hg pressure.

WARNING: To reduce the risks of leaking, slips and falls, and loss of therapy:

- Do not reuse. This product is designed for single patient use only.
- Do not use if end caps are not in place
- Ensure that all connections are securely tightened prior to priming set.

WARNING: To reduce the risk associated with patient exposure to leakage current that exceeds Type BF equipment limits:

- Do not use in combination with other medical equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment.

CAUTION: to reduce the risk of leaking and loss of therapy

- Do not remove warming cassette from warming unit while in use
- Do not exceed 245 mmHg pressure; ensure inline hand pressure pumps or other pumps do not exceed 245 mmHg pressure.

Figure 2

Instructions for Use

Do not use if package has been previously opened or damaged. Refer to the 3M™ Ranger™ irrigation fluid warming unit Model 247 Operator's Manual for set up and use of the warming system.

Set-up and Priming the irrigation disposable set

- Open warming set pouch.
- Ensure connectors are tight (Fig. 1, D, H)
- Present aseptically the inner pouch which contains the patient line (Fig. 1, I) to the sterile field.
- Before priming, slide the fluid warming cassette (Fig. 1, F) into the slot in the warming unit (Fig. 2)
- Close white pinch clamps (B1-B3, G1-2).
- Obtain the patient line (Fig. 1, I) and twist off end caps from irrigation set and patient line (Fig. 1, H).
- Connect the patient line to the irrigation set (Fig. 1, H).
- Insert the Ranger set spikes (Fig. 1, A) into the irrigation bag ports.
- To begin priming open all clamps. Ensure roller clamp (Fig. 1, J) is open.
- Invert the drip chamber (Fig. 1, C) until it fills at least 3/4 full, and then turn upright.
- Note:** Remove residual air bubbles by pinching the tubing.
- Close the white clamp (Fig. 1, G2) on the patient line when the line is free of air.
- Connect the patient line to the scope, Foley catheter, or similar device.
- Raise I.V. pole or pressurize fluid if indicated and open clamps when ready for use
- Turn the warming unit on.

The irrigation fluid warming set is now ready for use.

Removing the Irrigation warming set

- Turn the Ranger irrigation warming unit OFF.
- Close the blue inlet clamp (Fig. 1, E) proximal to the fluid warming cassette.
- Open all clamps distal to the cassette.
- Allow a small amount of fluid to drain from the end of the patient line to reduce pressure in the cassette.
- Remove the irrigation fluid warming cassette from the warming unit and discard according to institutional protocol (Fig. 2).

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Figure 1: Illustration of Disposable Irrigation Set

- | | |
|---------------------------------|----------------------------------|
| A. Spikes | G1-G2. White outlet pinch clamps |
| B1-B3. White inlet pinch clamps | H. Connector |
| C. Drip chamber | I. Patient Line |
| D. Connectors | J. Roller clamp |
| E. Blue inlet pinch clamp | K. Patient/scope connector |
| F. Fluid warming cassette | |

Symbol Glossary

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in Medical Device Regulation (EU) 2017/745 formerly EU Directive 93/42/EEC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. Source : ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Sterilized using ethylene oxide		Indicates a medical device that has been sterilized using ethylene oxide. Source: ISO 15223, 5.2.3
Do not use if package is damaged or open		Indicates a medical device that should not be used if the package has been damaged or opened. Source: ISO 15223, 5.2.8
Sterile fluid path		Indicates the presence of a sterile fluid path within the medical device in cases when other parts of the medical device, including the exterior, might not be supplied sterile. Source: ISO 15223, 5.2.9
Do not re-use		Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. Source: ISO 15223, 5.4.2
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Natural rubber latex is not present		Indicates natural rubber or dry natural rubber latex is not present as a material of construction within the medical device or the packaging of a medical device. Source: ISO 15223, 5.4.5 and Annex B
Medical Device		Indicates the item is a medical device.
Unique device identifier		Indicates bar code to scan product information into patient electronic health record
Importer		Indicates the entity importing the medical device into the EU.
Single sterile barrier system		To identify a single layer that makes up the sterile barrier system (SBS) and differentiate it from layers of protective packaging designed to prevent damage to the sterile barrier system and its contents

CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)

For more information see, [HCBGRegulatory.3M.com](https://www.hcbgregulatory.3m.com)

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3M and Ranger are trademarks of 3M. Used under license in Canada.
© 2025, 3M. Unauthorized use prohibited. All rights reserved.
3M Deutschland GmbH