
Ranger™

en	Irrigation Fluid Warming Unit, Model 247 Operator’s Manual	01
fr	Manuel de l'utilisateur de l'unité de réchauffement pour liquide d'irrigation, modèle 247	12
de	Erwärmungseinheit für Spülflüssigkeit, Modell 247, Bedienungsanleitung	24
it	Manuale d'uso dell'Unità riscaldante per fluidi di irrigazione, modello 247	36
es	Manual del operario de la unidad de calentamiento de fluidos de irrigación (Modelo 247)	48
nl	Bedieningshandleiding irrigatievloeistofverwarmingsunit, model 247	60
sv	Irrigationsvätska värmeenhet, modell 247 Operatörshandbok	71
da	Varmeenhed til irrigationsvæske, model 247, betjeningsmanual	82
no	Varmeenhet for irrigasjonsvæske, modell 247 Brukerhåndbok	93
fi	HuuhTELunesteen lämmitysyksikkö, malli 247, käyttöopas	104
pt	Manual do Operador da unidade de aquecimento de fluidos de irrigação, Modelo 247	114
el	Μονάδα Θέρμανσης Υγρού Καταιόνησης, Μοντέλο 247 Εγχειρίδιο Χειριστή	126
pl	Urządzenie do podgrzewania płynów stosowanych przy irygacji, model 247 — instrukcja obsługi	138
hu	Irrigációsolyadék-melegítő egység, 247-es modell, Kezelői kézikönyv	150
cs	Návod k obsluze ohřivací jednotky pro irigační tekutiny, model 247	161
sk	Používateľská príručka jednotky na ohrev irigačných tekutín, model 247	172
sl	Enota za ogrevanje tekočine za izpiranje, model 247 Uporabniški priročnik	183
et	Loputusvedeliku soojendussüsteemi mudeli 247 kasutusjuhend	194
lv	Skalošanas šķidrums sildītājs, modelis 247 – lietotāja rokasgrāmata	205
lt	Drėkinimo skysčio šildymo įrenginio (247 modelio) naudotojo vadovas	216
ro	Manualul operatorului unității de încălzire a fluidelor de irigare, model 247	227
ru	Инструкция по эксплуатации устройства для подогрева инфузионных растворов, модель 247	239
hr	Jedinica za zagrijavanje tekućine za ispiranje, model 247, priručnik za operatere	252
bg	Ръководство за оператора на устройство за затопляне на течности при иригация, модел 247	262
tr	İrigasyon Sıvısı Isıtma Ünitesi, Model 247 Kullanım Kılavuzu	274
zh	247 型输液加温器操作手册	284
ar	دليل المشغل لوحدة احتراز سوائل الإرواء، الطراز 247	293

 solventum
Internal Production Team

Solventum Production Artist:
Brian Pedersen
bnpedersen@solventum.com

Division:
MedSurg

This artwork has been produced by the
Solventum Internal Production Team

Artwork Stage:
7

Date:
9-Oct-2025

File Name: 34200292729_0264581_1.indd
Structure Spec # N/A

Artwork ID: 0264581.1
Material Spec #: N/A
Requestor: Debra Staack
PID #: 26060_s11

FOR PRINTER USE ONLY: This barcode symbol was produced using a bar width adjustment of -0.001". If this does not meet the needs of the printing process or equipment to be used, the printer is responsible to recreate the file using the appropriate bar width adjustment while meeting Solventum quality requirements.

Limitation of Liability: The supplied materials should be thoroughly inspected immediately upon receipt by the recipient for all content, placement and copy accuracy. It is also the responsibility of the recipient printer to examine and proof read all supplied materials before producing image carriers and printing final product.

Should the supplied materials be found to be incorrect please contact an Solventum representative immediately. The liability of Solventum shall be limited to correcting its own product. Solventum will not be responsible for any loss or damage including but not limited to the cost of product, print, image carriers, press delays, or other ancillary cost.

PROOF REQUIRED

Section 1: Technical Service and Order Placement..... 2
Technical Service and Order Placement..... 2
USA..... 2
Outside of the USA..... 2
In-warranty repair and exchange..... 2
When you call for technical support..... 2

Section 2: Introduction 2
Product Description 2
Indications for Use 2
Patient Population and Settings 2
Explanation of Signal Word Consequences 2
WARNING:..... 2
CAUTION:..... 3
NOTICE:..... 3
Overview and Operation..... 3
Ranger Irrigation Fluid Warming Set Example..... 3
Model 247 Product Safety Features..... 4

Section 3: Instructions for Use..... 4
Preparation and Setup of the Ranger Irrigation Fluid Warming Unit 4
Removing the Irrigation Fluid Warming Set from the Ranger Irrigation Fluid Warming Unit 4

Section 4: Troubleshooting 5

Section 5: General Maintenance and Storage 5
Cleaning Instructions 5
To Clean the Exterior of the Warming Unit:..... 5
To Clean the Heating Plates:..... 6
Method..... 6
To Clean Resistant, Dried-On Fluids:..... 6
Storage..... 6
Servicing 6
Symbol Glossary 6

Section 6: Specifications 7
Physical Characteristics..... 10
Classification..... 10
Environmental Conditions 10

Section 1: Technical Service and Order Placement

Technical Service and Order Placement

USA: TEL: 1-800-228-3957 (USA Only)

Outside of the USA: Contact your local 3M representative.

3M assumes no responsibility for the reliability, performance, or safety of the temperature warming unit if the following events occur:

- Modifications or repairs are not performed by a qualified, medical equipment service technician who is familiar with good practice for medical device repair.
- The unit is used in a manner other than that described in the Operator’s or Preventive Maintenance Manual.
- The unit is installed in an environment that does not provide grounded electrical outlets.
- The warming unit is not maintained in accordance with the procedures described in the Preventive Maintenance Manual.

In-warranty repair and exchange: To return a device to 3M for service, first obtain a return authorization (ra) number from a customer service representative. Please use the (ra) number on all correspondence when returning a device for service. A shipping carton will be delivered to you at no charge, if needed. Call your local supplier or sales representative to inquire about a loaner device while your device is being serviced.

When you call for technical support: We will need to know the serial number of your unit when you call us. The serial number label is located on the bottom of the warming unit.

Section 2: Introduction

Product Description

The Ranger irrigation fluid warming system consists of a Model 247 warming unit and a sterile disposable fluid warming set.

The 3M™ Ranger™ irrigation fluid warming system is designed to warm irrigation fluids and deliver these at flow rates under gravity up to 580 mL/min. Fluids can be warmed from 20°C inlet temperature to 33°C - 41°C. It takes less than 2 minutes to warm up to the 41°C set point temperature. An over-temperature alert point is set at 48°C and the under-temperature alert point is set at 33°C.

The Ranger irrigation fluid warming unit is designed to be mounted to an I.V. pole. A handle located on the top of the unit makes transport easy. When mounted to the I.V. pole, the unit fits easily above the 3M™ Bair Hugger™ warming unit. For additional information on Ranger disposable sets visit us online at rangerfluidwarming.com. This manual includes operating instructions and unit specifications for the Model 247 irrigation fluid warming system. For information about using Ranger irrigation fluid warming disposable sets with the Model 247 irrigation fluid warming unit, refer to the “Instructions for Use” included with each disposable component. The Ranger irrigation fluid warming system should only be used in healthcare facilities by trained medical professionals.

Indications for Use

The Ranger Irrigation Fluid Warming System is intended to warm irrigation fluids.

Patient Population and Settings

Adult and pediatric patients being treated in operating rooms, emergency trauma settings, or other areas where irrigation fluid warming is required.

Explanation of Signal Word Consequences

WARNING: Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.

NOTICE: Indicates a situation which, if not avoided, could result in property damage only.

WARNING:

1. To reduce the risks associated with hazardous voltage and fire and thermal energy hazards:
 - Do not substitute other devices (i.e. Model 245) for the Ranger irrigation fluid warming unit or Ranger irrigation disposable sets.
 - Do not continue use of the unit if the over-temperature alert continues to sound and the temperature does not return to the set point temperature. Immediately stop fluid flow and discard the disposable set. Have the warming unit tested by a biomedical technician or call 3M.
2. To reduce the risks associated with hazardous voltage and fire:
 - Do not modify or service this device, and do not open the warming unit case- as there are no user-serviceable parts in the Ranger irrigation fluid warming unit.
 - Connect power cord to receptacles marked “Hospital Only,” “Hospital Grade,” or a reliably grounded outlet.
 - Use only the 3M-supplied power cord intended for this product and certified for the country of use.
 - Do not allow the power cord to get wet.
 - Do not use the Ranger irrigation fluid warming system when it appears the warming unit, power cord, or warming set is damaged. Use only 3M specified replacement parts.
 - Keep power cord visible and accessible at all times. The plug on the power cord serves as the disconnect device. The wall socket outlet shall be as close as practical and shall be easily accessible.
 - Do not use with a multiple socket outlet or extension cord.
 - Do not use this equipment adjacent to, stacked with or in combination with other equipment without verifying that the total leakage current from the combined equipment does not exceed the safety limits for Type BF equipment and to assure normal operation in the configuration in which it will be used.
3. To reduce the risks associated with incorrect routing of fluids:
 - Ensure all luer connections are tightened.
4. To reduce the risk of injury associated with electromagnetic interference:
 - Care should be used when operating the Ranger irrigation fluid warming unit adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Ranger irrigation fluid warming unit and the other equipment should be observed to verify normal operation in the configuration in which it will be used.

CAUTION:

- 1. Not for IV use.
- 2. To reduce the risks associated with cross-contamination:
 - The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.
- 3. To reduce the risks associated with impact and facility medical device damage:
 - Clamp the Ranger irrigation fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm).
- 4. To reduce the risks associated with environmental contamination:
 - Follow applicable regulations when disposing of this device or any of its electronic components.
- 5. The Ranger irrigation fluid warming system has been tested to be resistant to both electromagnetic fields (EMI) and electrostatic discharge (ESD). To reduce the risk associated with EMI due to portable and mobile RF communications equipment:
 - Install and put into service the Ranger irrigation fluid warming system according to the EMC information provided in the Guidance and Manufacturer's Declaration.
 - Should interference occur move away from the portable or mobile RF communications equipment.
- 6. Under certain conditions, the Ranger irrigation fluid warming unit can become sensitive to ESD or surge, which can cause the audible alert to incorrectly activate. If such an event occurs, switch off the unit, then unplug and replug in the power cord from the power outlet to return to normal operation.

NOTICE:

- 1. The Ranger irrigation fluid warming unit meets medical electronic interference requirements. If radio frequency interference with other equipment should occur, it may be necessary to take mitigation measures such as reorienting, relocating, or connecting the warming unit to a different power source.
- 2. To avoid device damage:
 - Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not immerse the warming unit in cleaning or sterilizing solutions. The unit is not liquid proof.
 - Do not insert metallic instruments in the warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.

Overview and Operation

The Ranger irrigation fluid warming system consists of a Model 247 warming unit and a sterile disposable fluid warming set.

The warming unit is a compact, lightweight, liquid-resistant device with a clamp located on the side for attachment to an I.V. pole (see Figure 1). A carrying handle on the top of the unit makes it easy to transport.

Located on the front panel you will find:

- Alphanumeric display that indicates the heater temperature during normal operation. In an over temperature condition, the display alternately flashes a temperature of 48°C or higher and the word "HI." An audible alert also sounds. In an under-temperature condition, the display alternately flashes a temperature of 33°C or lower and the word "LO."
- Alert indicator light that comes on when either an over- or under-temperature condition occurs.


Figure 1
The Ranger irrigation fluid warming system consists of a model 247 warming unit and sterile warming set.

Ranger Irrigation Fluid Warming Set Example

Irrigation fluid warming sets¹ include a cassette, tubing, connectors, flow chamber and clamps (see Figure 2).

- A. Spikes
- B1-B5. White pinch clamps
- C. Drip chamber
- D. Connector
- E. Blue pinch clamp
- F. Fluid warming cassette
- G. Connector
- H. Patient line
- I. Roller clamp
- J. Scope connector


Figure 2

¹Refer to the instructions provided with each disposable set for information on use.

²With the bag hung 100 cm above the scope under gravity flow.

Model 247 Product Safety Features

The following chart describes the safety alarm features of the Ranger irrigation fluid warming unit.

Alert Type	Set-Point	Indication	Description
Primary over-temperature alert system	48°C	LED alternately displays "HI" and temperature, alert sounds, alert light flashes.	Primary Alert: Alert sounds, controller cuts power to the heaters until temperature falls below 48°C. When heater temperature drops below 48°C, audible and visual alerts stop.
Independent back-up over-temperature alert system (cut-off).	50°C	LED displays the temperature or has gone dark, alert sounds.	Primary failure. Temperature will rise to 50°C within a few seconds. The safety backup system is activated at 50°C and power is immediately shut off to the heating plates.
Under-temperature alert	33°C	LED displays "LO" and temperature, alert sounds, alert light flashes.	Heater temperature has fallen to 33°C. Alerts will stop when temperature rises above 33°C. Continue use of unit.

NOTE: The difference in over-temperature alert system cut-off points does not indicate any difference in safety and is simply a result of tighter controller tolerances.

The Ranger irrigation fluid warming system is designed with very tight controls to regulate its temperature to the 41°C set-point. The microprocessor controller measures the temperature at the fluid pathway inside the warming plate four times every second, making it very responsive to changes in flow rate.

Section 3: Instructions for Use

Preparation and Setup of the Ranger Irrigation Fluid Warming Unit

- 1. Attach the Ranger irrigation fluid warming unit to the I.V. pole. Tighten the pole clamp securely.

CAUTION:

To prevent tipping, clamp the Ranger irrigation fluid warming unit to an I.V. pole with a minimum 14" (35.6 cm) radius wheelbase and at a height no higher than 44" (112 cm). Failure to do so may result in damage to the product (see Figure 3).

- 2. Slide the irrigation cassette into the slot in the warming unit. The cassette can only fit into the device one way (see Figure 4, #1).
- 3. Prime the warming set. For more information about priming the set, refer to instructions provided with the warming sets (see Figure 4, #2).
- 4. Plug the power cord into an appropriate outlet. Turn the unit ON (see Figure 4, #3). In a few seconds the alpha-numeric display will illuminate. It takes less than 2 minutes to warm up to the 41°C set point temperature.
- 5. Begin flow. When flow is completed, remove the warming set and discard according to institutional protocol.

Removing the Irrigation Fluid Warming Set from the Ranger Irrigation Fluid Warming Unit

- 1. Close the inlet clamp proximal to the fluid warming cassette and open all clamps distal to the cassette.
- 2. Disconnect the warming set from the fluid source, if applicable.
- 3. Allow a small amount of fluid to drain from the end of the patient line to reduce pressure in the cassette.
- 4. Remove the warming set from the warming unit and discard according to institutional protocol.


Figure 3


Figure 4

Section 4: Troubleshooting

Condition	Cause	Solution
Nothing illuminates on the warming unit panel.	Unit is not turned on, plugged in, or power cord is not plugged into an appropriate outlet.	Turn unit on. Make sure the power cord is plugged into the power entry module of the warming unit. Make sure the warming unit is plugged into a properly grounded outlet.
	Unit failure.	Check panel fuses. Call 3M technical service.
Alert indicator illuminates and alert sounds, alpha-numeric display alternately flashes a temperature of 48°C or higher and the word “HI.”	Temporary over-temperature condition because:	
	An extreme change in flow rates occurred (e.g., from 870 ml/min to stop flow).	Open flow to reduce temperature. Alerts will stop when the temperature drops below 48°C. The unit is ready to use.
	Unit was turned on and reached set point temperature before warming cassette was inserted.	Alerts will stop when the temperature drops below 48°C. The unit is ready to use.
	Fluids were prewarmed to above 42°C before being run through the warming unit.	Turn off unit and unplug it. Discontinue infusion of fluids. Do not warm fluids before infusing them through the Ranger irrigation fluid warming unit.
Alert sounds, alpha-numeric display and alert indicator light go dark.	Primary controller failure. Unit will no longer operate.	Power to heating plates will shut off if temperature rises to 50°C. Turn unit off and unplug it. Discontinue use of unit. Discard disposable set. Alert will continue to sound if you do not unplug unit. Call 3M technical service.
Unit alerts soon after plugging it in (unit does not have to be turned on for this condition to occur). Heater temperature rises to 50°C and unit shuts down soon after plugging it in (unit does not have to be turned on for this condition to occur).	Test screw on bottom of unit is loose or missing.	Make sure test screw is completely tightened. If it is missing, turn unit off and unplug it. Call 3M technical service.
Alert sounds but warming unit has been turned off.	Independent backup safety system has been activated.	Unplug unit. Call 3M technical service.
Cannot remove cassette from unit.	Cassette is too full, fluids are still being infused, or clamp is open proximal to the warming cassette.	Make sure fluid is drained from warming cassette before sliding out the cassette, that fluids are no longer infusing, and that clamp is closed proximal to the warming cassette.
	Warming unit is below patient level, creating excessive back pressure.	Raise unit above patient level.
Alerts indicator light illuminates and alert sounds, alpha-numeric display alternately flashes a temperature of 33°C or below and the word “LO.”	Under temperature condition caused by very high flow using very cold fluid, or defective heater/relay.	Alerts should stop when temperature rises above 33°C. If alerts continue, turn off unit, unplug unit, and discontinue use. Call 3M technical service.
Alpha-numeric display reads “Er 4” or “Open.”	Open wire on temperature sensor.	Do not use unit. Call 3M technical service.
Alpha-numeric display reads “Er 5” or “Open.”	Electrical interference.	Remove the unit. Refer to biomedical technician or call 3M technical service.

Section 5: General Maintenance and Storage

Cleaning Instructions

To Clean the Exterior of the Warming Unit:

1. Disconnect the warming unit from the power source before cleaning.
2. Cleaning should be performed in accordance with hospital practices for cleaning OR equipment. After every use; wipe the warming unit and any other surfaces that may have been touched. Use a damp, soft cloth and a hospital approved mild detergent, germicidal disposable wipes, disinfecting towelettes, or antimicrobial spray. The following active ingredients are acceptable for use in cleaning the warming unit:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)
3. Let air dry.

CAUTION

To reduce the risks associated with cross-contamination:

- The cleaning tool provides only superficial cleaning, it does not disinfect or sterilize the interior of the unit.

To Clean the Heating Plates:

The Ranger hardware cleaning tool is intended to clean both heating plates of the warming unit. It is not necessary to disassemble the warming unit to use the tool.

Method

1.

Unplug the irrigation fluid warming unit.
2.

Unfold the cleaning tool. Wet the foam pads with a nonabrasive solution. The following active ingredients are acceptable for use with the cleaning tool:
 - Oxidizers (e.g. 10% Bleach)
 - Quaternary Ammonium Compounds (e.g. 3M™ Quat Disinfectant Cleaner)
 - Phenolics (e.g. 3M™ Phenolic Disinfectant Cleaner)
 - Alcohols (e.g. 70% Isopropyl Alcohol)
3.

Insert the tool from the back of the unit and pull the tool all the way out from the front (see Figure 5).
4.

Rinse the tool with water and repeat 3 times. Discard the tool according to institutional protocol.
5.

Wipe off the unit to remove excess fluid.


Figure 5

To Clean Resistant, Dried-On Fluids:

1.

Spray a nonabrasive solution inside the slot of the warming unit and let sit for 15-20 minutes.
2.

Clean the unit by using the cleaning tool.

NOTICE

1.

To avoid device damage:
 - Do not immerse the Ranger irrigation fluid warming unit or accessories in any liquid or subject them to any sterilization process.
 - Do not clean the warming unit with solvents. Damage to the case, label, and internal components may result.
 - Do not insert metallic instruments in the warming unit.
 - Do not use abrasive materials or solutions to clean the heater plates.
 - Do not allow spills to dry inside the unit, as this may make it more difficult to clean the unit.

NOTE: You may use a nonmetal instrument, such as a cotton swab, to clean the upper channels. If you are unable to adequately clean the unit, call 3M.

Storage

Store all components in a cool, dry place when not in use.

Servicing

There are no user-serviceable parts in the Ranger irrigation fluid warming unit. All service must be performed by 3M or an authorized service technician. In the USA, call 3M at 1-800-228-3957 (USA Only) for service information. Outside of the USA contact your local 3M representative.

Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Symbol Glossary

The following symbols may appear on the product’s labeling or exterior packaging.

Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer. Source: ISO 15223, 5.1.1
Authorized Representative in European Community / European Union		Indicates the authorized representative in the European Community / European Union. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Catalogue number		Indicates the manufacturer’s catalogue number so that the medical device can be identified. Source : ISO 15223, 5.1.6
Serial number		Indicates the manufacturer’s serial number so that a specific medical device can be identified. Source: ISO 15223, 5.1.7
Keep dry		Indicates a medical device that needs to be protected from moisture. Source: ISO 15223, 5.3.4
Caution		Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences. Source: ISO 15223, 5.4.4
Medical Device		Indicates the item is a medical device. Source: ISO 15223, 5.7.7
Unique device identifier		Indicates a carrier that contains unique device identifier information. Source: ISO 15223, 5.7.10
Importer		Indicates the entity importing the medical device into the locale. Source: ISO 15223, 5.1.8

CE Mark		Indicates conformity to all applicable European Union Regulations and Directives with notified body involvement.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a healthcare professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
UL Classified		Indicates product was evaluated and Listed by UL for the USA and Canada.
Refer to instruction manual/ booklet		To signify that the instruction manual/booklet must be read. ISO 7010-M002
"OFF" (power)		To indicate disconnection from the mains, at least for main switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5008
"ON" (power)		To indicate connection to the mains, at least for mains switches, or their positions, and all those cases where safety is involved. Source: IEC 60417-5007
Protective earth (ground)		To identify any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode. Source: IEC 60417, 5019
Equipotentiality		To identify the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential. Source: IEC 60417-5021
Fuse		Indicates a replaceable fuse
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle. Source: Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)
Type BF applied part		To identify a type BF applied part complying with IEC 60601-1. Source: IEC 60417-5333

For more information see, [HCBGregulatory.3M.com](https://www.hcbgregulatory.3m.com)

Section 6: Specifications

WARNING: Portable and mobile radio frequency (RF) communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Ranger Irrigation fluid warming unit, including cables specified by 3M. Otherwise, degradation of the performance of this equipment could result.

WARNING: The use of accessories and cables other than those specified or provided by 3M could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

EMC declaration tables according to IEC 60601-1-2:2014+A1:2020

Guidance and manufacturer's declaration – electromagnetic emissions

The model 247 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 247 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model 247 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer’s declaration – electromagnetic immunity

The model 247 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 247 should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV power lines	±2 kV power lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to gnd	±1 kV line to line ±2 kV line to gnd	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	70% U _T voltage dip for 30 cycles (@60 Hz) 0% U _T voltage dip for 1 cycle 0% U _T voltage dip for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0% U _T voltage interruption for 300 cycles (@60 Hz)	70% U _T voltage dip for 30 cycles (@60 Hz) 0% U _T voltage dip for 1 cycle 0% U _T voltage dip for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° phase angles 0% U _T voltage interruption for 300 cycles (@60 Hz)	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Immunity to proximity RF fields IEC 61000-4-3	See Radiated RF IEC 61000-4-3 Table below	See Radiated RF IEC 61000-4-3 Table below	See Radiated RF IEC 61000-4-3 Table below
Immunity to proximity magnetic fields IEC 61000-4-39	65 A/m, 134.2 kHz, 2.1 Khz pulse modulation	6.5 A/m	--
	7.5 A/m, 13.56 MHz, 50 kHz modulation	7.5 A/m	
Radiated RF electromagnetic fields and electromagnetic disturbances created by HF surgical equipment EC 60601-2-2:2017	Test guidance from IEC 60601-2-2:2017 Annex BB, Clause BB.4	Test guidance from IEC 60601-2-2:2017 Annex BB, Clause BB.4	HF surgical equipment operating at 100% power within proximity of the Ranger irrigation fluid warming unit and power cord.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer’s declaration – electromagnetic immunity

The model 247 is intended for use in the electromagnetic environment specified below. The customer or the user of the model 247 should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the model 247, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

- NOTE 1At 80 MHz and 800 MHz, the higher frequency range applies.
- NOTE 2These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model 247 is used exceeds the applicable RF compliance level above, the model 247 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model 247.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Compliance and test levels, Radiated RF IEC 61000-4-3

Test Frequency (MHz)	Band (MHz)	Modulation	Separation Distance (m)	Immunity Test Level (V/m)
385	380-390	Pulse, 18 Hz	0.3	27
450	430-470	FM ± 5 kHz, deviation 1 Khz sine	0.3	28
710	704-787	Pulse, 217 Hz	0.3	9
745		Pulse, 217 Hz	0.3	9
780		Pulse, 217 Hz	0.3	9
810	800-960	Pulse, 18 Hz	0.3	28
870		Pulse, 18 Hz	0.3	28
930		Pulse, 18 Hz	0.3	28
1720	1700-1990	Pulse, 217 Hz	0.3	28
1845		Pulse, 217 Hz	0.3	28
1970		Pulse, 217 Hz	0.3	28
2450	2400-2570	Pulse, 217 Hz	0.3	28
5240	5100-5800	Pulse, 217 Hz	0.3	9
5500		Pulse, 217 Hz	0.3	9
5785		Pulse, 217 Hz	0.3	9

Recommended separation distances between portable and mobile RF communications equipment and the model 247

The model 247 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the model 247 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the model 247 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.3 \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.30
10	3.69	3.69	7.37
100	11.67	11.67	23.30

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

- NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Physical Characteristics

Warming Unit

4.5 in. (11 cm) high x 7.5 in. (19 cm) wide x 10 in. (25 cm) long; wt.: 7 lb. 7 oz. (3.4 kg)

Classification

- Protection against electric shock: Class I Medical Electrical Equipment with Type BF applied part.
- Protection against ingress of water: IPX0 (Ordinary Equipment).
- Mode of operation: Continuous operation.
- Not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.

MEDICAL – GENERAL MEDICAL EQUIPMENT AS TO ELECTRIC SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH AAMI ES 60601-1:2005/(R)2012 + A1:2012/(R)2012 + A2:2021 and CAN/CSA-C22.2 No. 60601-1:14 (Reaffirmed 2022); Control No. 4HZ8

Electrical Characteristics	Temperature Characteristics and Accuracy
Input Voltage 100-120 or 220-240 VAC	Set Point Temperature 41 ± 1.5°C
Operating Frequency 100-120 VAC, 50/60 Hz 220-240 VAC, 50/60 Hz	Over-Temperature Alert 48 + 3/- 2°C
Maximum Heating Power 900 W	Under-Temperature Alert 33 ± 2°C
Fuse 2 x T10A-H (250V) for 100-120 VAC 2 x T6.3A-H (250V) for 220-240 VAC	Over-Temperature Cutoff 50 +2/-1°C

Fuse Type

Time delay, high breaking

Leakage Current

Meets leakage current requirements in accordance with IEC 60601-1.

Environmental Conditions

Operating Temperature Range

15° to 40°C (59° to 104°F)

Storage and Transport Temperature Range

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Section 1 : Assistance technique et commandes

Assistance technique et commandes

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3M, Ranger and Bair Hugger are trademarks of 3M.
All rights reserved.
3M Deutschland GmbH