Anaesthesia

(Coaxial Bain System) Mapleson D Modification

EN	Anaesthetic Circuit (Coaxial Bain System)	NO	Anestesikrets (Coaxial Bain System)
AR	دائرة التخدير (نظام Bain الم حوري)	PL	Układ anestetyczny (współosiowy układ Bain)
DE	Narkosegerät (Bain Koaxial-Beatmungssystem)	PT	Circuito anestésico (sistema de Bain coaxial)
ES	Circuito anestésico (sistema coaxial Bain)	PT-BR	Circuito anestésico (sistema Bain coaxial)
FI	Anestesian hengityspiiri (samankeskinen Bain-järjestelmä)	RO	Circuit anestezic (sistem coaxial Bain)
FR	Circuit d'anesthésie (système coaxial Bain)	RU	Анестезиологический контур (Коаксиальный контур Бэйна)
HU	Légzőkör (Koaxális csőrendszer)	SR	Anestetsko kolo (koaksijalni Bain sistem)
IT	Circuito per anestesia (Sistema coassiale Bain)	SV	Anestesikrets (koaxialt Bain-system)
JA	麻醉回路 (コアキシャル・ペインシステム) Mapleson D 改造	TH	วงจรสำหรับการให้ยาสลบ (ระบบ Bain แบบโคแอกเซียล) การปรับแต่ง Mapleson D
KO	마취 회로 (동축 베인 시스템) Mapleson D 개정	TL	Anaesthetic Circuit (Coaxial Bain System)
NL	Anesthesie circuit (coaxiaal Bainsysteem)	TR	Anestezi Devresi (Eş Eksenli Bain Sistemi)

EN	(English)	Anaesthetic Circuit (Coaxial Bain System) Mapleson D Modification	2
AR	(العربية)	دائرة التخدير (نظام Bain المحوري) Mapleson D تعديل	3
DE	(Deutsch)	Narkosegerät (Bain Koaxial-Beatmungssystem) Mapleson D	4
ES	(Español)	Circuito anestésico (sistema coaxial Bain) Modificación Mapleson D	5
FI	(Suomi)	Anestesian hengityspiiri (samankeskinen Bain-järjestelmä) Mapleson D -muunnos	6
FR	(Français)	Circuit d'anesthésie (système coaxial Bain) Modification Mapleson D	7
HU	(Magyar)	Légzőkör (Koaxális csőrendszer) Mapleson D módosítás	8
IT	(Italiano)	Circuito per anestesia (Sistema coassiale Bain) Modifica del sistema Mapleson D	9
JA	(日本語)	麻酔回路 (コアキシャル・ベインシステム) Mapleson D 改造	10
KO	(한국어)	마취 회로(동축 베인 시스템) Mapleson D 개정	11
NL	(Nederlands)	Anesthesie circuit (coaxiaal Bainsysteem) Mapleson D-modificatie	12
NO	(Norsk)	Anestesikrets (Coaxial Bain System) Mapleson D Modifikasjon	13
PL	(Język polski)	Układ anestetyczny (współosiowy układ Bain) Modyfikacja Mapleson D	14
PT	(Português)	Circuito anestésico (sistema de Bain coaxial) Modificação de Mapleson D	15
PT-BR	(Português Brasil)	Circuito anestésico (sistema Bain coaxial) Modificação Mapleson D	16
RO	(Română)	Circuit anestezic (sistem coaxial Bain) Modificare Mapleson D	17
RU	(Русский)	Анестезиологический контур (Коаксиальный контур Бэйна) модификации Mapleson D	18
SR	(Srpski)	Anestetsko kolo (koaksijalni Bain sistem) Mapleson D modifikacija	19
SV	(Svenska)	Anestesikrets (koaxialt Bain-system) Mapleson D-modifiering	20
TH	(ไทย)	วงจรสำหรับการให้ยาสลบ (ระบบ Bain แบบโคแอกเซียล) การปรับแต่ง Mapleson D	21
TL	(Tagalog)	Anaesthetic Circuit (Coaxial Bain System) Modipikasyon ng Mapleson D	22
TR	(Türkçe)	Anestezi Devresi (Eş Eksenli Bain Sistemi) Mapleson D Modifikasyonu	23

EN Anaesthetic Circuit (Coaxial Bain System) Mapleson D Modification

Key to Symbols:

Patient End

Adjust length if
required

Detachable
Connection

Instructions For Use

- This circuit is warranted for use on multiple anaesthesia patients during a period not exceeding 7-days from date of initial use. Use on subsequent patients in anaesthesia (to the first patient) must be in combination with a new (previously unused) suitably-validated bacterial viral filter/filter-HME located between the breathing system and the patient connection.
- Remove product from primary packaging and inspect contents for damage, ensuring that no separate components become mislaid.
- Complete any final assembly steps.
- Connect bacterial viral filter(s)/filter-HME (if appropriate).
- Assemble to anaesthesia machine (and to other system equipment, if required), ensuring that connections are secure and that unused ports are sealed.
- If a gas sampling line is connected, ensure that this is securely connected and monitored in use for occlusion by condensed water or by tube kinking.
- If an Adjustable Pressure Limiting (APL) valve is included with this circuit, use of the assembly is restricted to anaesthesia or anaesthesia reanimation only. For information on correct use of the APL valve and gas scavenging from the APL valve, please refer to specific IFU. Summary IFU for the APL valve is:
 - Rotate valve cap clockwise to apply increasing system pressure.
 - Rotate valve cap anti-clockwise to decrease system pressure and to facilitate spontaneously breathing patients.
 - Consider anaesthetic gas scavenging via the valve exhaust port.
- Verify flow patency and verify that separate inspiratory and expiratory gas paths are present and functioning.
- Perform other pre-use function test(s) as required for use on connected equipment.
- If the circuit does not pass these tests, verify correct set-up and re-test.
- In the event of repeated failure to pass the test(s), replace the system and/or contact a technician.
- When the system passes the test(s), enable therapy and connect to patient.
- Monitor performance of the system in use to ensure adequate and safe therapy.
- To avoid retention of CO₂, the following settings are recommended:
 - Spontaneous breathing:** set the APL valve to 'OPEN'. Adjust the total fresh gas flow to >120mL/Kg body weight per minute.
 - Controlled ventilation:** set the APL valve to partially-closed. Adjust the total fresh gas flow to >70mL/Kg body weight per minute.
- Ensure that use ceases within 7-days from date of initial use.

Warnings

- Keep package sealed until use – risk of blockage by small objects.
- A bacterial viral filter or filter-HME (if present) **must** be changed once in any 24-hour period or when it is suspected of interfering with effectiveness of ventilation.
- If reusing this system, it **must** be stored within a protective environment between uses. Alternatively, all ports and inlets, that are open to atmosphere, **must** be covered to protect them from exposure to contaminants.
- If using this circuit on a patient with known respiratory transmissible infection, dispose of after that use, regardless of air filtration strategy.
- Monitor for disconnections or other source of loss of effective ventilation therapy.
- Do not** wash, clean or attempt sterilisation.

Information

Device Name	Anaesthetic Circuit
Intended Use	Delivery of gaseous admixtures including oxygen and medical air from an anaesthesia machine to and from a patient
Indications for Use	Anaesthesia
Contraindications	Those patients not medically-indicated for the intended therapy
Warnings	Refer to any listed warnings
Intended Users	Appropriately trained healthcare professionals
Pre-use Precautions	Refer to Instructions For Use
Performance Characteristics	Refer to any listed performance characteristics
Disposal	Dispose as per the hospital's protocol for management of contaminated clinical waste
Device Malfunction	In the event of a device malfunction or an unexpected change in device performance, cease use of the device immediately
Serious Adverse Incidents	In the event of a malfunction of the device when in use on a patient and leading to a Serious Adverse Incident, contact the device manufacturer and the relevant National Competent Authority. (A Serious Adverse Incident is considered to be the death or serious deterioration in the state of health of the patient or user).
Use Restrictions	Use only with compatible devices. Consult Armstrong Medical for confirmation of compatible devices.
Other Information	This product may contain components affording antimicrobial protection (silver phosphate glass) to patients and users. See primary packaging for further information.

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