Resuscitation

Resuscitation set with APL valve

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EN Resuscitation Set with APL Valve (Mapleson Systems)

Key to Symbols:

Patient End

Adjust length if
required

Detachable
Connection

Instructions For Use

1. This product is warranted for use on 1 (one) intensive care patient or 1 anaesthesia patient for a maximum period of 7-days from date of initial use.
2. Remove product from primary packaging and inspect contents for damage, ensuring that no separate components become mislaid.
3. Complete any final assembly steps.
4. Connect bacterial viral breathing filter(s)/filter-HME (if appropriate).
5. Connect fresh gas tubing to gas supply source (oxygen flowmeter outlet or anaesthesia common gas outlet) and to other system equipment, if required, ensuring that connections are secure and that unused ports are sealed.
6. Prior to clinical use, occlude the patient connection with a cap or gloved thumb or forefinger. With fresh gas flow set at 10L/min and the APL valve cap rotated to the fully-closed position, observe that the reservoir bag inflates and distends for 3 (three) seconds. Then, rotate the APL valve cap anti-clockwise and observe that the reservoir bag collapses.
7. Verify flow patency alternately through the APL valve exhaust port and the patient connection.
8. Verify that separate inspiratory and expiratory gas paths are present and functioning.
9. Perform other pre-use function test(s) as required for use on connected equipment.
10. If the product does not pass these tests, verify correct set-up and re-test.
11. In the event of repeated failure to pass the test(s), replace the system and/or contact a technician.
12. When the product passes the test(s) connect to patient and commence therapy.
13. For information on correct use of the APL valve and gas scavenging from the APL valve, please refer to specific IFU. Summary IFU for the APL valve is:
 - a. Rotate valve cap clockwise to apply increasing system pressure.
 - b. Rotate valve cap anti-clockwise to decrease system pressure and to facilitate spontaneously breathing patients.
 - c. Consider anaesthetic gas scavenging via the valve exhaust port.
14. If a gas sampling line is connected, ensure that this is securely connected and monitored in use for occlusion by condensed water or by tube kinking.
15. Monitor performance of the system in use to ensure adequate and safe ventilation therapy.
16. Ensure that use ceases within 7-days from date of initial use.

Warnings

1. Keep package sealed until use – risk of blockage by small objects.
2. A bacterial viral breathing filter or filter-HME (if present) **must** be changed once in any 24-hour period or when it is suspected of interfering with effectiveness of ventilation.
3. Monitor for disconnections or other source of loss of effective ventilation therapy.
4. **Do not** wash, clean or attempt sterilisation.

Information

Device Name	Resuscitation Set (Mapleson Systems)
Intended Use	Delivery of oxygen from an oxygen flowmeter or anaesthesia machine common gas outlet to a patient for the purposes of anaesthesia, post-operative recovery, anaesthesia reanimation and critical care manual pulmonary hyperinflation
Indications for Use	Anaesthesia provision, anaesthesia reanimation or physiotherapy manual pulmonary hyperinflation
Contraindications	Those patients not medically-indicated for the intended therapy
Warnings	Refer to any listed warnings
Intended Users	Appropriately trained healthcare professionals
Pre-use Precautions	Refer to Instructions For Use
Performance Characteristics	Refer to any listed performance characteristics
Disposal	Dispose as per the hospital's protocol for management of contaminated clinical waste
Device Malfunction	In the event of a device malfunction or an unexpected change in device performance, cease use of the device immediately
Serious Adverse Incidents	In the event of a malfunction of the device when in use on a patient and leading to a Serious Adverse Incident, contact the device manufacturer and the relevant National Competent Authority. (A Serious Adverse Incident is considered to be the death or serious deterioration in the state of health of the patient or user).
Use Restrictions	Use only with compatible devices. Consult Armstrong Medical for confirmation of compatible devices.
Other Information	None declared.

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