Hi-Fatigue G Bone Cement

Instructions for use	2	en
Gebrauchsanweisung	8	de
Notice d'instruction	15	fr
Istruzioni per l'uso	22	it
Instrucciones de utilización	29	es
Instruções de utilização	36	pt
Gebruiksaanwijzing	43	nl
Οδηγίες χρήσης	50	el
Brugsanvisning	57	da
Bruksanvisning	63	sv
Kullanım talimatları	69	tr

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Instructions for use

Intended use and properties

Hi-Fatigue G Bone Cement is a gentamicin containing Poly(methyl methacrylate) based bone cement for fixation of endoprostheses in arthroplasty surgery. The added antibiotic, gentamicin sulfate, protects the implant and the surrounding tissue from colonization with pathogens that are sensitive to gentamicin.

Hi-Fatigue G Bone Cement consists of a powder and a liquid component that are mixed directly before use to form a ductile paste with optimal handling properties. After hardening Hi-Fatigue G Bone Cement provides a save and stable fixation of endoprostheses in orthopedic surgery. Hi-Fatigue G Bone Cement contains a radiopacifier to enable the identification of the cement in radiographs.

Composition Hi-Fatigue G Bone Cement

The cement powder contains

		1x20	1x40	1x60
Poly (methyl acrylate / methyl methacrylate)	27.54 %	5.75 g	11.5 g	17.25 g
Poly (methyl methacrylate / styrene)	57.49 %	12.0 g	24.0 g	36.0 g
Gentamicin (as gentamicin sulfate)	2.16 %	0.27 g (0.45 g)	0.55 g (0.9 g)	0.82 g (1.35 g)
Zirconium dioxide	11.98 %	2.5 g	5.0 g	7.5 g
Benzoyl peroxide	0.84 %	0.18 g	0.35 g	0.53 g

The liquid component contains

		1x20	1x40	1x60
Methyl methacrylate (stabilized with 60 ppm HQ)	99.35 %	9.93 ml	19.87 ml	29.8 ml
N,N-Dimethyl-p-toluidine	0.65 %	0.07 ml	0.13 ml	0.20 ml

Packaging sizes of Hi-Fatigue G Bone Cement

Pack size	Weight powder	Volume liquid
1x20	1 x 20.88 g	1 x 10 ml
2x20	2 x 20.88 g	2 x 10 ml
1x40	1 x 41.75 g	1 x 20 ml
2x40	2 x 41.75 g	2 x 20 ml
1x60	1 x 62.63 g	1 x 30 ml

Indications

Hi-Fatigue G Bone Cement is indicated for fixation of endoprostheses in arthroplasty surgery if an infection with organisms, that are sensitive to gentamicin, is diagnosed or suspected.

Contraindications

Hi-Fatigue G Bone Cement should not be used in patients with

- hypersensitivity to any of its constituents
- severe renal failure

Use during pregnancy and lactation

There are no adequate studies on the use of acrylate cements during pregnancy and lactation. The surgeon should weigh the benefit against the potential risk for mother and child before using Hi-Fatigue G Bone Cement.

Use in children

There are no adequate studies on the use of acrylate cements in children. The use of Hi-Fatigue G Bone Cement in patients who are still growing is inadvisable.

Side effects

In the course of every surgical or orthoplastic intervention, wound complications, superficial or deep wound infection, deep vein thrombosis as well as heterotopic bone formation and loosening or dislocation of the implant may occur.

After preparation of the implant site and immediately after application of the cement and insertion of the prosthesis, the rise in pressure in the medullary canal can lead to a temporary drop in blood pressure or short-term cardiac conduction irregularities. In rare cases, pulmonary embolism, myocardial infarction or apoplexy may occur. These cardiovascular and respiratory side effects which are known as implantation syndrome result mainly from infiltration of bone marrow constituents into the venous system.

In rare cases, the constituents of Hi-Fatigue G Bone Cement can cause allergic reactions or local irritations. There have been isolated reports of anaphylactic shock including cardiac arrest and sudden death.

Side effects specific of gentamicin are hearing problems and kidney damage. However, these side effects are extremely unlikely to occur because of the very low serum levels of gentamicin (< 1 µg/ml).

Monomer vapors can irritate the respiratory tract and eyes and may damage organs.

Interactions

No interactions of Hi-Fatigue G Bone Cement with other medical devices, drugs, materials, or techniques have been observed.

The tendency of gentamicin to block neuromuscular transmission can be intensified by concomitant administration of muscle relaxants, e.g. d-tubocurarine, suxamethonium or pancuronium, as well as by ether. Concomitant administration of potentially neurotoxic and/or nephrotoxic substances, e.g. cisplatin, other aminoglycosides, streptomycin, cefaloridine, viomycin, polymyxin B or polymyxin E, can increase the toxicity of gentamicin. However, interaction is very unlikely to occur because of the low serum levels of gentamicin.

Hi-Fatigue G Bone Cement does not pose a safety risk in a magnetic resonance environment.

Precautions

Instructions for handling, mixing and application of the cement and preparing the implant site must be followed carefully.

- Before using Hi-Fatigue G Bone Cement, the surgeon must be thoroughly familiar with its properties and handling characteristics. The handling and setting characteristics of Hi-Fatigue G Bone Cement depend on the temperature and mixing technique (see also graphics below). The surgeon is strongly recommended to carry out a trial run of the entire mixing, handling and setting process before performing a surgical procedure with Hi-Fatigue G Bone Cement. Inap-

appropriate handling of the cement may cause inadequate fixation leading to premature loosening or implant failure.

- If special mixing devices or application techniques are used, the surgeon must carefully read the relevant instructions first.
- Hardening of Hi-Fatigue G Bone Cement is an exothermic polymerization reaction. The released heat may damage bone or other tissues near to the cement.
- Methyl methacrylate is a slightly volatile, flammable liquid. The vapors produced during the mixing process can irritate the respiratory tract and eyes. Adequate ventilation and the use of closed mixing systems is recommended. Contact lenses should be protected from monomer vapors.
- Methyl methacrylate is lipid-soluble. Direct skin contact with the liquid monomer should be avoided. Rubber or latex gloves must not come into contact with the liquid component. It is advisable to wear an additional pair of polyethylene (PE) gloves under normal surgical gloves when handling the cement. The following materials have also proved to be suitable for protective gloves: PVP (polyethylene, ethylene vinyl alcohol, polyethylene) and Viton/Butyl.
- Gloves contaminated with cement mixture in the pre-dough state should be changed to avoid diffusion of monomer.
- Medical staff must inform the patient about all instructions and precautions to be taken. The time point of full weight bearing, as well as precautions and actions to be avoided after surgery should be judged by the responsible surgeon based on the surgical procedure, the type of associated implant and the condition of the patient.
- During and immediately after application of the cement and insertion of the prosthesis, blood pressure, pulse and respiration must be monitored carefully. If significant changes in cardiovascular or respiratory vital signs occur, appropriate measures are required.

Dosage

The required amount of Hi-Fatigue G Bone Cement depends on the anatomical conditions of the patient and the used prostheses. If a large amount is needed, up to four 1x40 units can be mixed together. The whole content of all units must be mixed. It is advisable to keep additional packs of Hi-Fatigue G Bone Cement ready as a reserve.

Preparation

Depending on the mixing technique, appropriate tools should be available. Sterile porcelain bowls, stainless steel bowls or plastic bowls that are suitable for monomers, sterile mixing spoons or spatulas made of porcelain or high-grade steel, or a sterile mixing system for bone cements.

Prechilling of the cement is not necessary.

The cement powder is packed in three nested pouches. The outer pouch and the outer side of the middle pouch are non-sterile. The inner side of the middle pouch and the inner pouch are sterile. The ampoule with monomer liquid is packed in a blister. The outer side of the blister is non-sterile while the inner side and the ampoule are sterile. The outer pouch and the blister should be removed from the carton in the non-sterile area. The outer and the middle pouch ("peel-off") as well as the blister should be opened by a nurse. The sterile inner pouch and the ampoule should be aseptically handed over to a member of the surgical team in the sterile area. The ampoule should be opened by breaking its neck using the attached protective cap. The inner pouch should be opened with sterile scissors.

Note to preparation


1. The powder and liquid components have been carefully designed to complement each other. The entire contents of the pouch and ampoule must always be mixed.

USING ONLY PARTS OF THE COMPONENTS IS NOT PERMITTED.

The addition of any other powders or liquids can reduce solidity and impair handling characteristics.

2. For preparation the temperature of the components must not be below 4 °C.
3. Temperature-Time Correlation: Hi-Fatigue G Bone Cement is designed for handling between 18 °C and 23 °C. The phases of mixing, waiting and application are correlated to the polymerization process which is strongly dependent on the temperature of the components and environment. Viscosity increases with the progression of polymerization. Setting, waiting and handling time is reduced by higher temperatures and prolonged by lower temperatures (see Temperature-Time Correlation graphic below).

Figure: Temperature-Time Correlation of HI-Fatigue G Bone Cement (not pre-chilled)

A: Mixing phase

B: Waiting phase

C: Application phase

D: Hardening phase

Vacuum mixing

The use of a vacuum mixing system is recommended, representing the state-of-the-art mixing technology which reduces air entrapment. The detailed instructions for mixing and application techniques with special devices can be found in the instructions for use supplied by the manufacturer of the mixing system.

The contents of the required ampoules are poured into the mixing vessel and the same number of powder units are added. All ampoules and pouches must be emptied thoroughly. The cement should be mixed for about 30 seconds.

Manual mixing

For manual mixing, the entire contents of the required ampoules are poured into a suitable sterile, inert mixing vessel. It is not recommended to prechill the bone cement before manual mixing. The entire contents of the same number of powder units are added to the liquid. The two components are carefully mixed using a suitable spatula until a homogenous paste has formed (this needs usually about 30 seconds). It is important to avoid entrapment of air during mixing and application.

Application of the cement

After mixing the cement should be rested for 1-2 minutes (Waiting phase, see the figure of Temperature-Time Correlation) which enables the cement to reach the appropriate viscosity. The end of waiting phase is indicated when the cement does not stick to the gloves any more. For application a cement syringe or other application systems can be used. The detailed instructions for using special application systems can be found in the instructions for use supplied by the manufacturer of the application system.

During the cement application and the following prosthesis insertion, the cardiovascular and respiratory function of the patient must be carefully monitored.

The application phase ends when the cement becomes rubbery and elastic, and no longer binds completely when kneaded. Cement at this stadium must not be applied any more. Otherwise, an optimal filling cannot be ensured and there is a risk of premature implant loosening.

Surplus cement should be removed before it hardens.

For an optimal fixation of the implant, it should be inserted within the application phase and held still in place until the cement hardens (Hardening phase, see the figure of Temperature-Time Correlation).

To avoid undesirable effects and to ensure long-term stable fixation of the implant, State-of-the-art cementing techniques should be applied. This includes careful preparation of the implant site with pulsatile lavage and drying before cement application, adequate drainage to prevent over pressure in the medullary canal during implantation, and the use of a cement restrictor to enhance cement penetration and prevent leakage. Over pressurizing of the bone cement is to be avoided as this may lead to extrusion of the bone cement beyond the site of the intended application and may damage the surrounding tissue. The entire medullary canal should be completely filled with cement to form a cement mantle with an ideal thickness of 2-5 mm. It is recommended to verify the correct implantation using suitable imaging procedures.

Removal of bone cement in case of revision

For revision, a radical debridement is necessary to remove the loosened prosthesis, the cement and devitalized tissues. Reaming of the medullary canal followed by a thorough lavation using a pulse lavage system is recommended. The procedure to remove the implant is a case by case decision by the responsible surgeon depending on the condition of the implant and the patient.

Storage

Hi-Fatigue G Bone Cement should be stored away from direct sunlight and at temperatures between 0 °C (32 °F) and +25 °C (77 °F). Do not use Hi-Fatigue G Bone Cement which has been exposed to temperatures exceeding the limitations during storage and transport.

Sterilization

The monomer is sterile filtered and filled into the ampoules under aseptic conditions. The blistered ampoules and the cement powder are sterilized with ethylene oxide.

The unused contents of any opened or damaged packs must not be re-sterilized and should be discarded.

Shelf life

The expiry date is printed on the batch label of the carton, the outer and inner pouches, the ampoule and the blister lid. Hi-Fatigue G Bone Cement must not be used after the expiry date. The polymer powder must not be used if it shows yellowish discoloration.

Single-use

Hi-Fatigue G Bone Cement must never be reused. Due to the functionality (hardening) of PMMA bone cements, Hi-Fatigue G Bone Cement is only suitable for use within the specified application period. Once opened, any pack should be completely used or discarded and not retained for use on a later occasion. One unit of Hi-Fatigue G Bone Cement must be used for only one single patient.

Safe disposal

1. Mixed cement should be discarded with hospital waste only after fully hardening.
2. For separate disposal of powder and liquid, please ask local waste disposal authorities.

Information

For further information, please contact your supplier or the manufacturer directly.

Symbol definitions

	Manufacturer		Do not re-use
	Use-by date		Consult instructions for use
	Batch code		Caution
	Catalogue number		CE marking of conformity
	Sterilized using aseptic processing techniques		Flammable
	Sterilized using ethylene oxide		Harmful
	Do not re-sterilize		Health Hazard
	Do not use if package is damaged		Ampoule
	Keep away from sunlight		Blister
	Temperature limit		Pouch

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