

Dear Valued Customer,

As part of Mölnlycke's ongoing commitment to providing the most effective products to healthcare professionals, we are transitioning the original style of Mepilex® Border and Mepilex® Border Lite to the upgraded Mepilex® Border Comfort range. Mepilex Border Comfort is an improved product range, of dressings designed specifically to stay in place and uniquely conform thanks to the superior combination of Flex Technology and Safetac®. By upgrading you will receive a familiar product with enhanced features for no additional costs across the range.

Mepilex® Border Comfort (5-layered bordered silicone foam dressing)

- **360° proprietary* Flex Technology allows Mepilex® Border Comfort to move in every direction, reducing skin stress, increasing comfort and keeping it in place^{1,2}.**
- **Smart Exudate Management:**
 1. Rapid absorption with minimal swelling: The foam layer allows rapid absorption of both normal and viscous exudate^{3,4}.
 2. **Improved spreading layer: Maximises capacity by channelling exudate. Exudate is channeled away from the wound and over a wider area to increase moisture vapor loss and reduce swelling^{6,7}.**
 3. **Traps bacteria away from the wound: The increased superabsorbent fibers in the retention layer prevent exudate-containing bacteria from re-entering the wound bed, even under compression⁸.** The backing film prevents the entry of external bacteria⁹.
 4. **Tracks exudate without dressing removal: The Exudate Progress Monitor allows the spread of exudate to be objectively assessed and recorded, without disturbing the wound, potentially avoiding excess dressing changes¹⁰.**
- Safetac® technology: for gentle dressing changes¹¹.

Mepilex® Border Comfort Lite (4-layered bordered silicone foam dressing)

- **Original 4 layers with the addition of Flex technology: 360° proprietary* Flex Technology allows Mepilex® Border Comfort Lite to move in every direction, reducing skin stress, increasing comfort and keeping it in place^{1,2}. Offering up to 77% better conformability than Mepilex® Border and Mepilex® Border Lite^{12,13}.**
- **Addition of the Exudate Progress Monitor allows Health Care Practitioners objectively track and record exudates, potentially avoiding excess dressing changes¹⁰.**
- Safetac® technology: for gentle dressing changes¹¹.

We appreciate your partnership as a valued Mölnlycke customer, and we would like to thank you for your continued business. For more information please contact your local representative or call Mölnlycke Health Care Customer Service T: 0800 917 4919 or email CSC.uk@Molnlycke.co.uk. ***NB: Improvements Bolded**

Thank you for your support.

Natalie Milne – UK Product Manager

References

Mölnlycke Health Care. Mepilex® Border Flex – External data – FEM simulations Mx8 Flex vs AllevynLife. Report no. PD-529747. Data on file. 2. Mölnlycke Health Care. FEM –Mepilex® Border Flex vs Aquacel foam. Report no. PD-534571. Data on file. 3. Mölnlycke Health Care. Mepilex® Border Flex – Absorption and dispersion tests on inclined plane, viscous test solution. Report no. PD-528871. Data on file. 4. Mölnlycke Health Care. Product manual Initial abs different test solutions. Report no. PD-535352. Data on file. 6. Mölnlycke Health Care. Mepilex® Border Flex Product Manual – Spreading layer PD-528874. Data on file. 7. Mölnlycke Health Care. Product Manual – Maintained dressing profile. Report no. 525458. Data on file. 8. Mölnlycke Health Care. Mepilex® Border Flex-Bacteria encapsulation. Report no. PD-537072. Data on file. 9. Mölnlycke Health Care. Product Manual Mepilex® Border Flex – Viral penetration. Report no. PD-535090. Data on file. 10. Mölnlycke Health Care. Mepilex® Border Flex – Estimate of spreading area using dot pattern on backing film. Report no. PD-528872. Data on file. 11. White, R. A multinational survey of the assessment of pain when removing dressings. Wounds UK 2008;4(1):14-22. 12. Mölnlycke Health Care. Data on file (2019). 13. Mölnlycke Health Care. Data on file (2019)