

URGENT: FIELD SAFETY NOTICE

BLUselect® Tracheostomy Tube Kits , BLUselect® Suctionaid® Tracheostomy Tube Kits, BLUgriggs® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with BLUselect® Tracheostomy Tube with or without Forceps, BLUperc® Dilation Procedural Tray with Single Stage Dilator Products, BLUperc® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with or without BLUselect® Tracheostomy Tube


15th August 2024

Dear Valued BLUSelect® Customers,

Smiths Medical, Inc. is issuing this Urgent Field Safety Notice to notify you of a potential defect with the following BLUSelect®, BLUgriggs® and BLUperc® products listed on pages 5-18. This letter details the issue and the required steps for you to complete.

Issue:

Smiths Medical has identified the potential for a disconnection of the pilot balloon from the tracheostomy inflation line within specific lots of the BLUSelect®, BLUgriggs® and BLUperc® products because of a manufacturing defect. See the photo example of the issue below.



Potential Risk:

If the pilot balloon used to inflate the tracheostomy cuff becomes detached from the inflation line, the cuff pressure may not be maintained which can lead to inadequate ventilation of, and increased risk of aspiration to the patient. To date, Smiths Medical has received thirteen (13) reports of serious injury and zero (0) deaths associated with this issue.

Affected Product

Please see the affected item and lot numbers on pages 5-18 and also product distribution date range.

Smiths Medical Actions:

Smiths Medical is sending this notification to all BLUSelect®, BLUgriggs® and BLUperc® customers who received products from Smiths Medical listed on pages 5-18. Smiths Medical will provide credit to affected customers upon receipt of a completed response form to certify product destruction.

Customer Required Actions:

When using the device, all instructions, including warning and cautions contained in the Instructions for Use Documentation must be followed with heightened awareness. Please complete the following actions listed below

1. Check all inventory locations within your institution for the affected catalog numbers and lot numbers listed in the notification and discontinue use. Discard all affected products following your institution's process for discarding. If discarding is not immediately possible at your facility, then the product should be quarantined until disposal.
2. Share this notification with all potential users of the device to ensure they are aware of this notification and proposed mitigations. If the devices are used at another location, please ensure this communication is delivered there.
3. Complete and return the attached Customer Response Form to EMEA-FSN@icumed.com within 10 days of receipt to acknowledge your understanding of this notification.
4. **DISTRIBUTORS:** If you have distributed potentially affected products to your customers, please immediately forward this notice to them and request that they complete the response form and return it to **YOU**. Then the **DISTRIBUTOR** must complete a SINGLE form with the required details and return to EMEA-FSN@icumed.com

For further inquiries, please contact Smiths Medical using the following information:

Smiths Medical Contact	Contact Information	Areas of Support
Global Complaint Management	globalcomplaints@icumed.com	To report adverse events or product complaints
Customer Service	https://www.icumed.com/about-us/contact-us	Questions about product replacement and/or credit.

Your country regulatory agency has been notified of this action

Smiths Medical is committed to patient safety and is focused on providing exceptional product reliability and the highest level of customer satisfaction. Thank you for your prompt support on this important matter. We appreciate your cooperation.

Sincerely,



Andy Mathein
Vice President of Quality

Enclosures:

- Customer Response Form (See Below)
- Affected Product List (see pages 5-18)

URGENT: FIELD SAFETY NOTICE – RESPONSE FORM

BLUselect® Tracheostomy Tube Kits , BLUselect® Suctionaid® Tracheostomy Tube Kits, BLUgriggs® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with BLUselect® Tracheostomy Tube with or without Forceps, BLUperc® Dilation Procedural Tray with Single Stage Dilator Products, BLUperc® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with or without BLUselect® Tracheostomy Tube

15th August 2024

Check your inventory and complete the information below, even if you do not have the affected product. *Failure to complete all sections of this page may result in improper, delayed or denied credit.*

Please return the completed form to EMEA-FSN@icumed.com, If you have questions about this form please contact EMEA-FSN@icumed.com or your local sales representative.

Name of Hospital / Facility	
Hospital / Facility Address	
Telephone Number	
Name and Title of Person Completing this Form	
Signature of Person Completing this Form	
Date	
If Purchased through a distributor, please list distributor name/location here for traceability purposes	

Please select one:

- I have **NO** affected products (complete and return this form to the e-mail address above)
- YES**, I have affected products, I have notified users in my facility and I have followed the instructions provided to me and destroyed all affected items (see table below)

If you have affected product on hand, please complete table 1 below:

TABLE 1

Item / SKU Number	Lot Number	Quantity in inventory (Eaches)	Quantity Destroyed (Eaches)	Date of Destruction

If you have distributed the product further, please complete table 2 below with collated information received from your customers and respond to ICU Medical with the overall information.

TABLE 2

Item / SKU Number	Lot Number	Quantity destroyed locally (Eaches)	Date of Destruction

Adverse events and complaints associated with the use of this product should be reported and emailed to Smiths Medical's Global Complaint Management Department at globalcomplaints@icumed.com.

Affected items below were distributed in the UK between Jan 2020 and Apr 2024

SKU/List/Model Number	Product Description	Lot Number
101/541/070	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 7.0mm BLUselect® trach tube	4202486
101/541/080	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 8.0mm BLUselect® trach tube	3912758
101/541/080	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 8.0mm BLUselect® trach tube	4170885
101/543/070	BLUgriggs® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Trach Tube, no forceps	3887541
101/543/070	BLUgriggs® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Trach Tube, no forceps	4310304
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	3887555
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	3928959
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	4005346
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	4217858
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	4234862
101/561/070	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 7.0 mm BLUselect trach tube	4348810
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	3887536
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	3906566
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	3906567
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	3906568
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4217859
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4232677
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4238252
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4281149
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4281150
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4283700
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4304152
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4348884

SKU/List/Model Number	Product Description	Lot Number
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4348885
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4348892
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4348893
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4382036
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4382040
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4386516
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4386518
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4396545
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4399796
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4399798
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4405137
101/561/080	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 8.0 mm BLUselect trach tube	4412103
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	3887528
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	3976736
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4170968
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4205592
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4234863
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4242949
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4281078
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4285761
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4310484
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4310487
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4310488
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4382034

SKU/List/Model Number	Product Description	Lot Number
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4397559
101/561/090	BLUperc Percutaneous Dilation Procedural Kit with Single Stage Dilator and 9.0 mm BLUselect trach tube	4397560
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	3925814
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4010058
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4017238
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4054722
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4054723
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4054801
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4165746
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4234758
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4245279
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4249227
101/563/070	BLUperc® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube	4281148
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	3910412
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	3925815
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4054726
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4070503
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4070504
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4100287
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4100288
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4103806
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4115434
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4147040
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4149715

SKU/List/Model Number	Product Description	Lot Number
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4175503
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4199278
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4205593
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4208994
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4217863
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4224557
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4234759
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4281042
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4283716
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4283717
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4283718
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4285934
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4285935
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348839
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348840
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348845
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348847
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348848
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348849
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348850
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348851
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4348860
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4350945
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4354160

SKU/List/Model Number	Product Description	Lot Number
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4385041
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4385043
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4385046
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4385050
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4396414
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4396420
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4396535
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4396536
101/563/080	BLUperc® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube	4396537
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	3980266
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4009974
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4061111
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4094349
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4115435
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4156962
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4170970
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4212065
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4227249
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4238178
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4281031
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4348865
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4348866
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4348867
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4386559

SKU/List/Model Number	Product Description	Lot Number
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4408641
101/563/090	BLUperc® Percutaneous Dilation Procedural Kit with 9.0mm BLUselect® Suctionaid® Trach Tube	4412148
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	3925835
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4025600
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4106716
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4230149
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4252500
101/800/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4382002
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	3895101
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4022790
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4037975
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4152440
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4156971
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4217804
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4238237
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4249341
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4305944
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4305947
101/800/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4368641
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4025599
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4037974
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4117927
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4152439
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4156970
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4238234
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4249234
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4348898
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4348899
101/800/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4365892
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	3895103
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	3933044
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4037971
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4052621
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4138225
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4156969
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4224695
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4230078
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4232686

SKU/List/Model Number	Product Description	Lot Number
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4245469
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4247407
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4270931
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4305933
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4305934
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4305935
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4305936
101/800/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4305938
101/800/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	3951221
101/800/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4034824
101/800/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4256383
101/800/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4348861
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	3951220
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4197025
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4217803
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4238236
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4261648
101/800/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4306047
101/800/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4230076
101/800/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4306058
101/800/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4372161
101/802/060CZ	BLUSELECT 6.0, CUFFED, FEN	3925839
101/802/060CZ	BLUSELECT 6.0, CUFFED, FEN	4028452
101/802/060CZ	BLUSELECT 6.0, CUFFED, FEN	4382046
101/802/070CZ	BLUSELECT 7.0, CUFFED, FEN	3925740
101/802/070CZ	BLUSELECT 7.0, CUFFED, FEN	3972531
101/802/075CZ	BLUSELECT 7.5, CUFFED, FEN	3968265
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	3951185
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	4268672
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	4310441
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	4423828
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	4436040
101/802/080CZ	BLUSELECT 8.0, CUFFED, FEN	4456953
101/802/085CZ	BLUSELECT 8.5, CUFFED, FEN	4028454
101/802/090CZ	BLUSELECT 9.0, CUFFED, FEN	3954847
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	3887495
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	3933190
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	3948313
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	3958315

SKU/List/Model Number	Product Description	Lot Number
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4106722
101/810/060CZ	BLUSELECT 6.0, CUFFED, NON-FEN	4268673
101/810/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	3887478
101/810/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4106813
101/810/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4243097
101/810/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4254571
101/810/070CZ	BLUSELECT 7.0, CUFFED, NON-FEN	4285856
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	3887504
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	3925795
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4010055
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4073763
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4177908
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4222083
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4249233
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4281153
101/810/075CZ	BLUSELECT 7.5, CUFFED, NON-FEN	4381994
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	3887481
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4173265
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4207825
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4217114
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4230143
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4285773
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4306016
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4312272
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4381980
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4387365
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4403461
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4415817
101/810/080CZ	BLUSELECT 8.0, CUFFED, NON-FEN	4423832
101/810/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4067915
101/810/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4152435
101/810/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4178022
101/810/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4281079
101/810/085CZ	BLUSELECT 8.5, CUFFED, NON-FEN	4306011
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	3993623
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4106736
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4173266
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4243117
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4306027

SKU/List/Model Number	Product Description	Lot Number
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4327669
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4403367
101/810/090CZ	BLUSELECT 9.0, CUFFED, NON-FEN	4458561
101/810/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4175497
101/810/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4217875
101/810/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4249333
101/810/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4310443
101/810/100CZ	BLUSELECT 10.0, CUFFED, NON-FEN	4387368
101/812/060CZ	BLUSELECT 6.0, CUFFED, FEN	3887479
101/812/060CZ	BLUSELECT 6.0, CUFFED, FEN	3959080
101/812/060CZ	BLUSELECT 6.0, CUFFED, FEN	3976752
101/812/060CZ	BLUSELECT 6.0, CUFFED, FEN	4106752
101/812/060CZ	BLUSELECT 6.0, CUFFED, FEN	4294444
101/812/070CZ	BLUSELECT 7.0, CUFFED, FEN	3887499
101/812/070CZ	BLUSELECT 7.0, CUFFED, FEN	3925781
101/812/070CZ	BLUSELECT 7.0, CUFFED, FEN	3968279
101/812/070CZ	BLUSELECT 7.0, CUFFED, FEN	4106754
101/812/070CZ	BLUSELECT 7.0, CUFFED, FEN	4242910
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	3887552
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	3968280
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	4217831
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	4268677
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	4381979
101/812/075CZ	BLUSELECT 7.5, CUFFED, FEN	4413594
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	3887488
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	3948295
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	3989255
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	4073770
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	4106756
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	4227246
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	4234757
101/812/080CZ	BLUSELECT 8.0, CUFFED, FEN	4385065
101/812/085CZ	BLUSELECT 8.5, CUFFED, FEN	4222085
101/812/085CZ	BLUSELECT 8.5, CUFFED, FEN	4261649
101/812/085CZ	BLUSELECT 8.5, CUFFED, FEN	4382043
101/812/090CZ	BLUSELECT 9.0, CUFFED, FEN	3972311
101/812/090CZ	BLUSELECT 9.0, CUFFED, FEN	4170983
101/812/090CZ	BLUSELECT 9.0, CUFFED, FEN	4217120
101/812/090CZ	BLUSELECT 9.0, CUFFED, FEN	4243118

SKU/List/Model Number	Product Description	Lot Number
101/812/100CZ	BLUSELECT 10.0, CUFFED, FEN	4064797
101/812/100CZ	BLUSELECT 10.0, CUFFED, FEN	4259283
101/812/100CZ	BLUSELECT 10.0, CUFFED, FEN	4413618
101/812/100CZ	BLUSELECT 10.0, CUFFED, FEN	4423838
101/812/100CZ	BLUSELECT 10.0, CUFFED, FEN	4443290
101/860/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4058593
101/860/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4117932
101/860/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4242923
101/860/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4376593
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	3887509
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	3959094
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4062929
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4067465
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4170834
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4186645
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4198020
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4202488
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4232682
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4234708
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4234709
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4252464
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4294329
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4303383
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4306007
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4306008
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4368707
101/860/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4370061
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4065369
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4189388
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4199270
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4224578
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4245286
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4285843
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4317594
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4386514
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4387362
101/860/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4408613
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	3895095
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4017244

SKU/List/Model Number	Product Description	Lot Number
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4054804
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4061062
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4071181
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4178020
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4191167
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4195032
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4202463
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4211996
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4211997
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4219600
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4232674
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4234707
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4242937
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4252471
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4283735
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4305958
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4312278
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4357891
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4384987
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4386575
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4388728
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4392809
101/860/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4405274
101/860/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4205591
101/860/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4214079
101/860/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4249237
101/860/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4270940
101/860/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4310532
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	3895091
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4037924
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4058605
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4191506
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4209002
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4239848
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4252466
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4285839
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4306064
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4310533
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4382009

SKU/List/Model Number	Product Description	Lot Number
101/860/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4385029
101/860/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4054789
101/860/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4199355
101/860/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4412135
101/870/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	3925734
101/870/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	3951306
101/870/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4106745
101/870/060CZ	BLUSELECT 6.0, SUCTIONAID, CUFFED, NON-FEN	4217882
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	3968404
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4022711
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4034792
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4058662
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4073762
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4138214
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4156953
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4165744
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4178001
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4217826
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4249291
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4256336
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4283719
101/870/070CZ	BLUSELECT 7.0, SUCTIONAID, CUFFED, NON-FEN	4291171
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4058606
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4061022
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4138176
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4183954
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4212004
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4238185
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4249236
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4270941
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4276631
101/870/075CZ	BLUSELECT 7.5, SUCTIONAID, CUFFED, NON-FEN	4382010
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	3887515
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4051928
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4054796
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4070522
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4117931
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4134866
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4138056

SKU/List/Model Number	Product Description	Lot Number
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4163106
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4208997
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4238168
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4254564
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4303380
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4310523
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4312285
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4317581
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4319800
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4368642
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4388736
101/870/080CZ	BLUSELECT 8.0, SUCTIONAID, CUFFED, NON-FEN	4390747
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	3887484
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4002986
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4054842
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4067915
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4085217
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4183981
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4212070
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4245285
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4249238
101/870/085CZ	BLUSELECT 8.5, SUCTIONAID, CUFFED, NON-FEN	4376585
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4000188
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4061185
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4106832
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4155508
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4178002
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4183974
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4242948
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4249307
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4249334
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4276656
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4285763
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4285764
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4303381
101/870/090CZ	BLUSELECT 9.0, SUCTIONAID, CUFFED, NON-FEN	4387367
101/870/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4075798
101/870/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4106718
101/870/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4209008

SKU/List/Model Number	Product Description	Lot Number
101/870/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4217806
101/870/100CZ	BLUSELECT 10.0, SUCTIONAID, CUFFED, NON-FEN	4306010
101/891/070	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 7.0mm BLUselect® Suctionaid® Trach Tube	4232678
101/891/080	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 8.0mm BLUselect® Suctionaid® Trach Tube	3887523
101/891/090	BLUgriggs® Percutaneous Dilation Procedural Kit with forceps and 9.0mm BLUselect® Suctionaid® Trach Tube	4310412
101/893/070	BLUgriggs® Percutaneous Dilation Procedural Kit with 7.0mm BLUselect® Suctionaid® Trach Tube, no forceps	4189353
101/893/080	BLUgriggs® Percutaneous Dilation Procedural Kit with 8.0mm BLUselect® Suctionaid® Trach Tube, no forceps	4310419