AURUM

This chart shows the needle-free connectors compatible with the Aurum Legacy and Aurum CONNECT ranges of pre-filled syringes.

Intended for healthcare professionals only. This is non-promotional material. This chart has been prepared and funded by ethypharm.

The following list of needleless luer access device connectors are compatible with the Aurum *Legacy* range of pre-filled syringes:

1 x 10ml Prefilled Syringe Adrenaline (Epinephrine) Injection 1:10,000



Baxter Healthcare - Clearlink and V-Link Connectors B.Braun Medical - Ultrasite/Safeflow/Caresite Connectors

1mg in 10ml

Becton Dickinson (BD) - Q-Syte Connectors

BMDi TUTA Healthcare's Luer valve Connectors, including: H-Port, N-Pulse and T-Site

BMS Critical Care - EasiFlush Connectors

CareFusion - Maxplus Displacement Connectors and SmartSite Valve Connectors

CODAN Ltd - SWAN-LOCK Connectors

Fannin - Flowart Connectors

GBUK Healthcare - Prosum NF001 Connectors

ICU Medical UK - Nuitiv Connectors and Tego D1000 Connectors

RPG Medical - EZ-IVNF0001 Connectors

Spirit Medical - SecureConnect Connectors

Vygon - Vadsite Connectors

The following needless luer access device connectors are compatible with

the Aurum range of pre-filled syrings when used with an adaptor:

Vygon - Bionector Connectors: with ATP431 adaptor. Contact Vygon (UK) Ltd for information on how to obtain the adaptors.


ICU Medical UK - MicroCLAVE and MicroCLAVE Clear and Neutron Connectors: with the use of the 011-CS25 Clear Adaptor. Contact ICU Medical UK for information on how to obtain the adaptors.

To differentiate between both barrels, there are clear visual differences in the packaging as follows:

- Text: Updated layout with the word "CONNECT."
- Icon: "Twist to open" icon on CONNECT packaging.
- Colour: "Green" with white stripes box closure label on CONNECT packaging (previously white).

The following list of needleless luer access device connectors are compatible with the Aurum CONNECT range of pre-filled syringes, without the use of an adaptor:

Adrenaline (Epinephrine) Injection 1:10,000 1 x 10ml Prefilled Syringe CONNECT For IV. IO or ET tube us 1mg in 10m



Baxter Healthcare - Clearlink/V-Link/OneLink Connectors B.Braun Medical - Ultrasite/Safeflow/Rymed Invision Plus/ Caresite Connectors

Becton Dickinson (BD) - Q-Syte/MaxZero/SmartSite/NeutraClear Connectors BMDi TUTA Healthcare's Luer valve Connectors, including: H-Port, N-Pulse and T-Site BMS Critical Care - EasiFlush Connectors CareFusion - Maxplus Displacement Connectors and SmartSite Valve Connectors **CODAN Ltd - SWAN-LOCK Connectors** Fannin - Flowart Connectors **GBUK Healthcare - Prosum NF001 Connectors** ICU Medical UK - Nuitiv/Tego D1000/MicroCLAVE/MicroCLAVE Clear/NanoCLAVE Clear/ **Neutron Connectors** Promedicare - B-site Connectors **RPG Medical - EZ-IVNF0001 Connectors** Sendal UK - TPS Connectors Spirit Medical - SecureConnect Connectors Vygon - Bionector & Vadsite Connectors

Quality Department Ref Number: DR24-285 (Version 1) Quality Department Ref Number: DR22-321 (Version 7)

Other connectors are not currently recommended