

Choosing the right blanket for your patient.



In order to ensure that your facility has access to the right forced-air warming blanket for your patients, Solventum has a wide range of options available to suit different positions and procedures.

What to remember when choosing a forced-air warming blanket:

- When using forced-air warming blankets or gowns to help maintain your patient's core temperature, the area covered by the blanket is of utmost importance. The larger the area covered, the more efficient the forced-air warming blanket will be.¹
- Forced-air warming blankets do not only transfer heat to the patient. They also reduce heat loss from the covered area to zero, which means the skin under the blanket is no-longer a source of heat loss but one of heat gain, changing the heat balance of the body.¹
- Ensure that when using all blankets, the perforated side is against the patient's skin.

Bair Hugger™ forced-air warming blankets and alternatives

The Bair Hugger™ Universal Warming Gown made with Thinsulate™ Insulation can be used as an alternative for any blanket and can be utilised across the perioperative pathway – Pre, Intra and Post. Sizes available: Small 81105, Medium 81106, Large 81107.

Adult full body					
Model 300 full body blanket	315	570	635	545	
Model 315 multi access full body blanket	300	570	635		
Adult lower body					
Model 525 lower body blanket	622	542	537	635	
Adult blanket					
Model 622 multi-position upper body 300	300	525	542	545	635
Model 542 dual port torso blanket	622	525	545	635	
Adult underbody series					
Model 545 adult full access underbody	635	585	550	622	542
Model 575 spinal underbody	525	570	622		
Model 585 lithotomy blanket	622	542	550		
Model 635 underbody blanket	585	545	542	550	622
Sterile cardiac					
Model 630 sterile cardiac blanket non-sterile alternatives	635	585	545	550	
Paediatric over body					
Model 537 small lower body blanket	525	622	550		
Paediatric underbody					
Model 550 paediatric large full access underbody blanket	585	545	575	635	542
Model 555 paediatric small full access underbody blanket	550				
Specialty and complex surgery					
Model 570 surgical access blanket	525	622	542	635	545

Request a consultation with a Solventum representative or Clinical Specialist to discuss blanket options for your patients at solventum.com/en-gb/home/medical/surgical-solutions/temperature-management/