

Information for Clinical Choice

A review of Vaginal Pessaries

Support Document

Introduction

Information for Clinical Choice (ICC) has been developed to assist Clinicians in the decision-making process when assessing the suitability of a product by providing a clear illustration and description of the features of a range of similar products supplied through NHS Supply Chain.

The Supporting Documents are the result of a product review, conducted by the Clinical Engagement and Implementation Managers (CEIM), with support from clinical stakeholders across the NHS, as part of our assurance process.

The aim, alongside delivering savings back into NHS frontline services, is to ensure that clinical choice remains at the forefront of any product switching decision.

Maternity, Obstetrics & Gynaecology: Pessaries

Vaginal pessaries are a non-surgical (conservative) method of managing pelvic organ prolapse and are used intravaginally to try to restore the prolapsed organs to their normal position and relieve symptoms¹. Choice and fit of pessary are based on clinician experience, and the fitting process is described as 'trial and error' until a type and size are found that is comfortable for the woman during all activities¹.

There is currently a diverse portfolio of pessaries available via NHS Supply Chain to meet a variety of clinical indications, therefore the tool below has been developed to help Clinicians and Buyers compare ranges based on product type, size, and indications for use. In addition to this, product descriptions have been updated to identify those products containing metal as this has safety implications for patients undergoing MRI investigations.

Products have also been price ranked to facilitate a simpler price comparison.

Reference List

1 Pelvic, Obstetric & Gynaecological Physiotherapy (POGP) & United Kingdom Continence Society (UKCS) (2021) *UK Clinical Guideline for best practice in the use of vaginal pessaries for pelvic organ prolapse*
(<https://www.ukcs.uk.net/resources/Documents/Pessary%202021/UK%20Pessary%20Guideline%20final%20April21.pdf>)



If you have any questions, would like further information, or have feedback to share, please contact:




Jo Norris – jo.norris@supplychain.nhs.uk



Clinical Engagement & Implementation Manager, Maternity, Obstetrics & Gynaecology



Claire Addy – claire.addy@supplychain.nhs.uk



Category Specialist, Maternity, Obstetrics & Gynaecology



Cerclage / Modified Cup		
Supplier	RB Medical	RB Medical
NPC	Size 1 65 32 21mm (FVC1777) Size 2 70 32 21mm (FVC1739) Size 3 65 35 21mm (FVC1643) Size 4 70 35 21mm (FVC1757) Size 5 65 35 30mm (FVC1705) Size 6 70 35 30mm (FVC1737)	Size 1 65 32 21mm (FVC1663) Size 2 70 32 21mm (FVC1651) Size 3 65 35 21mm (FVC1702) Size 4 70 35 21mm (FVC1802) Size 5 65 35 30mm (FVC1645) Size 6 70 35 30mm (FVC1807)
Description	Cerclage / Modified Cup with drain	Cerclage / Modified Cup without drain
Picture		
UOI	1	1
Route	eDirect	eDirect
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Not included	Not included
Indications for Use (Manufacturer's)	Mild Pelvic Organ Prolapse	Mild Pelvic Organ Prolapse
Country of Manufacture	USA	USA




Cube			
Supplier	RB Medical	Mediplus	Mediplus
NPC	Size 0 25mm (FVC1786) Size 1 29mm (FVC1712) Size 2 32mm (FVC1774) Size 3 36mm (FVC1700) Size 4 40mm (FVC1801) Size 5 45mm (FVC1713) Size 6 49mm (FVC1684) Size 7 55mm (FVC1753) Size 8 63mm (FVC1616) Size 9 69mm (FVC1791) Size 10 74mm (FVC1680)	Size 0 25mm (FVC1286) Size 1 30mm (FVC241) Size 2 35mm (FVC242) Size 3 38mm (FVC243) Size 4 41mm (FVC245) Size 5 44mm (FVC246) Size 6 51mm (FVC1889) Size 7 57mm (FVC1113)	38x30x68mm (FVC1892) 41x35x76mm (FVC1186) 44x38x82mm (FVC1187) 51x44x92mm (FVC1898) 57x44x101mm (FVC1188)
Description	Cube	Cube	Cube Tandem
Picture			
UOI	1	1	1
Route	eDirect	Blue Diamond	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Severe Pelvic Organ Prolapse	Severe Pelvic Organ Prolapse	Severe Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	3 rd degree Pelvic Organ Prolapse Procidentia Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Stress Urinary Incontinence Cystocele Mild Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Stress Urinary Incontinence Cystocele Mild Rectocele
Country of Manufacture	USA	USA	USA



Cube with drain		
Supplier	RB Medical	Mediplus
NPC	Size 0 25mm(FVC1693) Size 1 29mm (FVC1736) Size 2 32mm (FVC1697) Size 3 36mm (FVC1787) Size 4 40mm (FVC1652) Size 5 45mm (FVC1808) Size 6 49mm (FVC1785) Size 7 55mm (FVC1635) Size 8 63mm (FVC1769) Size 9 69mm (FVC1754) Size 10 74mm (FVC1707)	Size 0 25mm (FVC1114) Size 1 30mm (FVC1155) Size 2 35mm (FVC1156) Size 3 38mm (FVC1099) Size 4 41mm (FVC1157) Size 5 44mm (FVC1158) Size 6 51mm (FVC1159) Size 7 57mm (FVC1160)
Description	Cube with drain	Cube with drain
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Severe Pelvic Organ Prolapse	Severe Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	3 rd degree Pelvic Organ Prolapse Procidentia Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Stress Urinary Incontinence Cystocele Mild Rectocele
Country of Manufacture	USA	USA



Cup		
Supplier	RB Medical	RB Medical
NPC	Size 0 48mm (FVC1629) Size 1 52mm (FVC1677) Size 2 57mm (FVC1799) Size 3 63mm (FVC1612) Size 4 68mm (FVC1657) Size 5 73mm (FVC1740) Size 6 77mm (FVC1622) Size 7 82mm (FVC1708)	Size 0 48mm (FVC1678) Size 1 52mm (FVC1806) Size 2 57mm (FVC1756) Size 3 63mm (FVC1778) Size 4 68mm (FVC1634) Size 5 73mm (FVC1719) Size 6 77mm (FVC1695) Size 7 82mm (FVC1607)
Description	Cup	Cup with support
Picture		
UOI	1	1
Route	eDirect	eDirect
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Not included	Not included
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse	1 st -2 nd degree Pelvic Organ Prolapse
Country of Manufacture	USA	USA

Dish		
Supplier	RB Medical	Mediplus
NPC	Size 0 50mm (FVC1667) Size 1 55mm (FVC1654) Size 2 60mm (FVC1696) Size 3 63mm (FVC1779) Size 4 68mm (FVC1775) Size 5 73mm (FVC1648) Size 6 80mm (FVC1803) Size 7 85mm (FVC1764)	Size 0 55mm (FVC286) Size 1 60mm (FVC289) Size 2 65mm (FVC290) Size 3 70mm (FVC291) Size 4 75mm (FVC292) Size 5 80mm (FVC1273) Size 6 85mm (FVC294) (all sizes contain metal)
Description	Dish	Dish
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Stress Urinary Incontinence Large anterior Pelvic Organ Prolapse +/- Stress Urinary Incontinence	Stress Urinary Incontinence Large anterior Pelvic Organ Prolapse +/- Stress Urinary Incontinence
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele Urinary Incontinence	1 st -2 nd degree Pelvic Organ Prolapse Cystocele Stress Urinary Incontinence
Country of Manufacture	USA	USA



Dish with support		
Supplier	RB Medical	Mediplus
NPC	Size 0 50mm (FVC1762) Size 1 55mm (FVC1810) Size 2 60mm (FVC1617) Size 3 63mm (FVC1742) Size 4 68mm (FVC1703) Size 5 73mm (FVC1748) Size 6 80mm (FVC1664) Size 7 85mm (FVC1609)	Size 0 55mm (FVC1274) Size 1 60mm (FVC306) Size 2 65mm (FVC1897) Size 3 70mm (FVC308) Size 4 75mm (FVC309) Size 5 80mm (FVC1884) Size 6 85mm (FVC1275) (all sizes contain metal)
Description	Dish with support	Dish with support
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Stress Urinary Incontinence Large anterior Pelvic Organ Prolapse +/- Stress Urinary Incontinence	Stress Urinary Incontinence Large anterior Pelvic Organ Prolapse +/- Stress Urinary Incontinence
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele Urinary Incontinence	1 st -2 nd degree Pelvic Organ Prolapse Mild Cystocele Stress Urinary Incontinence
Country of Manufacture	USA	USA





Donut			
Supplier	Medicare Colgate	RB Medical	Mediplus
NPC	Size 3 64mm (FVC1716) Size 4 70mm (FVC1794) Size 5 76mm (FVC1596) Size 6 83mm (FVC1639)	Size 0 51mm (FVC1641) Size 1 57mm (FVC1682) Size 2 64mm (FVC1679) Size 3 70mm (FVC1766) Size 4 76mm (FVC1718) Size 5 83mm (FVC1597) Size 6 89mm (FVC1666) Size 7 95mm (FVC1610)	Size 0 51mm (FVC271) Size 1 57mm (FVC272) Size 2 64mm (FVC274) Size 370mm (FVC1161) Size 4 76mm (FVC1595) Size 5 83mm (FVC279) Size 6 89mm (FVC280) Size 7 95mm (FVC285)
Description	Donut	Donut	Donut
Picture			
UOI	1	1	1
Route	Blue Diamond	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Severe Pelvic Organ Prolapse	Severe Pelvic Organ Prolapse	Severe Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	Pelvic Organ Prolapse Stress Urinary Incontinence	2 nd -3 rd degree Pelvic Organ Prolapse Cystocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Mild Cystocele Rectocele
Country of Manufacture	Sweden	USA	USA


Gehrung		
Supplier	Mediplus	Mediplus
NPC	Size 0 27x38x35mm (FVC1162) Size 1 29x44x38mm (FVC1219) Size 2 35x51x41mm (FVC1899) Size 3 38x54x44mm (FVC232) Size 4 41x57x48mm (FVC1164) Size 5 44x60x51mm (FVC1165) Size 6 48x64x54mm (FVC1166) Size 7 51x73x57mm (FVC1167) Size 8 54x76x60mm (FVC1100) Size 9 57x79x64mm (FVC1168) (all sizes contain metal except Size 0)	Size 0 27x38x35mm (FVC1115) Size 1 29x44x38mm (FVC1116) Size 2 35x51x41mm (FVC1117) Size 3 38x54x44mm (FVC1118) Size 4 41x57x48mm (FVC458) Size 5 44x60x51mm (FVC1890) Size 6 48x64x54mm (FVC1903) Size 7 51x73x57mm (FVC1894) Size 8 54x76x60mm (FVC1119) Size 9 57x79x64mm (FVC1120) (all sizes contain metal except size 0)
Description	Gehrung	Gehrung with knob
Picture		
UOI	1	1
Route	Blue Diamond	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	2 nd -3 rd degree Pelvic Organ Prolapse Stress Urinary Incontinence Cystocele Mild Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Stress Urinary Incontinence Cystocele Mild Rectocele
Country of Manufacture	USA	USA



Gellhorn short stem with drain		
Supplier	RB Medical	Mediplus
NPC	Size 1 44mm (FVC1650) Size 2 51mm (FVC1615) Size 3 57mm (FVC1706) Size 4 63mm (FVC1795) Size 5 70mm (FVC1676) Size 6 76mm (FVC1717) Size 7 82mm (FVC1611) Size 8 89mm (FVC1636)	38mm (FVC764) 44mm (FVC1577) 51mm (FVC1124) 57mm (FVC767) 64mm (FVC1585) 70mm (FVC1586) 76mm (FVC1581) 83mm (FVC1583) 89mm (FVC1582) 95mm (FVC1588)
Description	Gellhorn short stem with drain	Gellhorn short stem with drain
Picture		
UOI	1	1
Route	eDirect	eDirect
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	More advanced Pelvic Organ Prolapse	More advanced Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	2 nd -3 rd degree Pelvic Organ Prolapse / Procidentia Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Mild Cystocele Rectocele
Country of Manufacture	USA	USA



Gellhorn standard stem with drain			
Supplier	RB Medical	Mediplus	Mediplus
NPC	Size 1 44mm (FVC1728) Size 2 51mm (FVC1763) Size 3 57mm (FVC1733) Size 4 63mm (FVC1683) Size 5 70mm (FVC1776) Size 6 76mm (FVC1689) Size 7 82mm (FVC1796) Size 8 89mm (FVC1598)	38mm (FVC145) 44mm (FVC153) 51mm (FVC1121) 57mm (FVC160) 64mm (FVC161) 70mm (FVC168) 76mm (FVC1122) 83mm (FVC170) 89mm (FVC171) 95mm (FVC172)	51mm (FVC1191) 57mm (FVC259) 64mm (FVC265) 70mm (FVC266) 76mm (FVC1192) 83mm (FVC268) 89mm (FVC269)
Description	Gellhorn standard stem with drain	Gellhorn standard stem with drain	Gellhorn standard stem with drain
Picture			
UOI	1	1	1
Route	eDirect	Blue Diamond	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Rigid Silicone
Indications for Use ¹	More advanced Pelvic Organ Prolapse	More advanced Pelvic Organ Prolapse	More advanced Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	2 nd -3 rd degree Pelvic Organ Prolapse / Proctentia Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Mild Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Vaginal Vault Mild Cystocele Rectocele
Country of Manufacture	USA	USA	USA


Hodge		
Supplier	RB Medical	Mediplus
NPC	Size 0 65mm (FVC1722) Size 1 70mm (FVC1747) Size 2 75mm (FVC1685) Size 3 80mm (FVC1602) Size 4 85mm (FVC1731) Size 5 90mm (FVC1726) Size 6 95mm (FVC1653) Size 7 100mm (FVC1649) Size 8 105mm (FVC1752) Size 9 110mm (FVC1804) (all sizes contain metal)	Size 4 92x54mm (FVC1887) (contains metal)
Description	Hodge	Hodge
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele Stress Urinary Incontinence Small Vaginal Introitus	Stress Urinary Incontinence Mild Cystocele Cervical Incompetence Retrodisplacement
Country of Manufacture	USA	USA




Hodge				
Supplier	RB Medical	Mediplus	Mediplus	Mediplus
NPC	Size 0 65mm (FVC1691) Size 1 70mm (FVC1783) Size 2 75mm (FVC1723) Size 3 80mm (FVC1655) Size 4 85mm (FVC1800) Size 5 90mm (FVC1788) Size 6 95mm (FVC1690) Size 7 100mm (FVC1751) Size 8 105mm (FVC1704) Size 9 110mm (FVC1688) (all sizes contain metal)	Size 0 70x44mm (FVC1153) Size 2 83x49mm(FVC1893) Size 3 89x51mm (FVC1902) Size 5 95x57mm (FVC1908) Size 6 98x60mm (FVC1900) Size 7 108x64mm (FVC1895) Size 8 117x67mm (FVC1891) (all sizes contain metal)	Size 0 70mm x 44mm (FVC1133) (contains metal)	Size 0 70mm x 44mm (FVC1143) (contains metal)
Description	Hodge with support	Hodge with support	Hodge with knob	Hodge with knob & support
Picture				
UOI	1	1	1	1
Route	eDirect	Blue Diamond	Blue Diamond	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele Stress Urinary Incontinence Exercise Incontinence Small Vaginal Introitus	1 st -2 nd degree Pelvic Organ Prolapse Stress Urinary Incontinence Mild Cystocele	Stress Urinary Incontinence	Stress Urinary Incontinence
Country of Manufacture	USA	USA	USA	USA




Inflatable Donut	
Supplier	Mediplus
NPC	51mm (FVC127) 57mm (FVC140) 64mm (FVC141) 70mm (FVC144)
Description	Inflatable Donut
Picture	
UOI	1
Route	Blue Diamond
Material	Latex
Indications for Use ¹	Severe Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	2 nd -3 rd degree Pelvic Organ Prolapse Mild Cystocele Mild Rectocele
Country of Manufacture	USA




Marland		
Supplier	RB Medical	RB Medical
NPC	Size 2 57mm (FVC1675) Size 3 63mm (FVC1699) Size 4 70mm (FVC1698) Size 5 76mm (FVC1797) Size 6 83mm (FVC1681) Size 7 89mm (FVC1805) Size 8 95mm (FVC1749)	Size 2 57mm (FVC1781) Size 3 63mm (FVC1771) Size 4 70mm (FVC1673) Size 5 76mm (FVC1789) Size 6 83mm (FVC1686) Size 7 89mm (FVC1773) Size 8 95mm (FVC1660)
Description	Marland	Marland with support
Picture		
UOI	1	1
Route	eDirect	eDirect
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	2 nd -3 rd degree Pelvic Organ Prolapse Moderate to Extensive Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Moderate to Extensive Cystocele Rectocele
Country of Manufacture	USA	USA




Oval		
Supplier	RB Medical	RB Medical
NPC	Size 1 51mm (FVC1626) Size 2 57mm (FVC1671) Size 3 64mm (FVC1658) Size 4 69mm (FVC1735) Size 5 76mm (FVC1674) Size 6 83mm (FVC1765) Size 7 89mm (FVC1603) Size 8 95mm (FVC1750) Size 9 102mm (FVC1687)	Size 1 51mm (FVC1792) Size 2 57mm (FVC1637) Size 3 64mm (FVC1710) Size 4 69mm (FVC1647) Size 5 76mm (FVC1720) Size 6 83mm (FVC1759) Size 7 89mm (FVC1665) Size 8 95mm (FVC1631) Size 9 102mm (FVC1632)
Description	Oval	Oval with support
Picture		
UOI	1	1
Route	eDirect	eDirect
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Not included	Not included
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse +/- Cystocele (designed to fit narrow vaginal vault)	2 nd -3 rd degree Pelvic Organ Prolapse Moderate to Extensive Cystocele Rectocele
Country of Manufacture	USA	USA



Regula	
Supplier	Mediplus
NPC	Size 0 25x49x35mm (FVC1276) Size 1 27x51x38mm (FVC1277) Size 2 29x54x41mm (FVC1885) Size 3 32x56x44mm (FVC1905) Size 4 35x57x48mm (FVC1896) Size 5 37x60x51mm (FVC1278) Size 6 38x67x54mm (FVC1279) Size 7 41x73x57mm (FVC1280) Size 8 44x79x60mm (FVC1281) (all sizes contain metal except Size 0 and 1)
Description	Regula
Picture	
UOI	1
Route	Blue Diamond
Material	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse
Country of Manufacture	USA


Ring			
Supplier	Bray Group	Cooper Surgical	Bray Group
NPC	50mm (FVC1835)	50mm (FVC483)	50mm (FVC1853)
	53mm (FVC1831)	53mm FVC1591)	53mm (FVC1857)
	56mm (FVC1817)	56mm (FVC485)	56mm (FVC1839)
	59mm (FVC1848)	59mm (FVC486)	59mm (FVC1844)
	62mm (FVC1814)	62mm (FVC487)	62mm (FVC1846)
	65mm (FVC1830)	65mm (FVC488)	65mm (FVC1833)
	68mm (FVC1843)	68mm (FVC489)	68mm (FVC1592)
	71mm (FVC1813)	71mm (FVC490)	71mm (FVC1852)
	74mm (FVC1816)	74mm (FVC491)	74mm (FVC1858)
	77mm (FVC1829)	77mm (FVC492)	77mm (FVC1847)
	80mm (FVC1819)	80mm (FVC493)	80mm (FVC1812)
	85mm (FVC1856)	85mm (FVC494)	85mm (FVC1821)
	90mm (FVC1820)	90mm (FVC495)	90mm (FVC1832)
	95mm (FVC1818)	95mm (FVC496)	95mm (FVC1834)
	100mm (FVC1855)	100mm (FVC497)	100mm (FVC1854)
	110mm (FVC1845)	110mm (FVC498)	
	Description	Ring	Ring
Picture			
UOI	20	20	10
Route	Blue Diamond	eDirect	Blue Diamond
Material	Flexible PVC	Flexible Vinyl	Rigid Polythene
Indications for Use ¹	Any type of Pelvic Organ Prolapse	Any type of Pelvic Organ Prolapse	Any type of Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	Pelvic Organ Prolapse	Pelvic Organ Prolapse	Pelvic Organ Prolapse
Country of Manufacture	UK	USA	UK



Ring			
Supplier	Medicare Colgate	RB Medical	Mediplus
NPC	70mm (FVC503)	Size 1 51mm (FVC1627) Size 2 57mm (FVC1780) Size 3 62mm (FVC1725) Size 4 70mm (FVC1661) Size 5 74mm (FVC1811) Size 6 82mm (FVC1613) Size 7 89mm (FVC1600) Size 8 96mm (FVC1694) Size 9 102mm (FVC1638)	Size 0 44mm (FVC1282) Size 1 51mm (FVC1283) Size 2 57mm (FVC048) Size 3 64mm (FVC049) Size 4 70mm (FVC1594) Size 5 76mm (FVC051) Size 6 83mm (FVC052) Size 7 89mm (FVC055) Size 8 95mm (FVC056) Size 9 102mm (FVC057) Size 10 108mm (FVC058) Size 11 114mm (FVC060) Size 12 121mm (FVC1886) Size 13 127mm (FVC1284) (sizes 11, 12 and 13 contain metal)
Description	Ring	Ring	Ring
Picture			
UOI	1	1	1
Route	Blue Diamond	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Any type of Pelvic Organ Prolapse	Any type of Pelvic Organ Prolapse	Any type of Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele	Mild 1 st degree Pelvic Organ Prolapse	1 st -2 nd degree Pelvic Organ Prolapse Cystocele
Country of Manufacture	China	USA	USA



Ring with support			
Supplier	Medicare Colgate	RB Medical	Mediplus
NPC	76mm (FVC677)	Size 1 51mm (FVC1721) Size 2 57mm (FVC1744) Size 3 62mm (FVC1790) Size 4 70mm (FVC1761) Size 5 74mm (FVC1784) Size 6 82mm (FVC1709) Size 7 89mm (FVC1646) Size 8 96mm (FVC1770) Size 9 102mm (FVC1601)	Size 0 44mm (FVC082) Size 1 51mm (FVC085) Size 2 57mm (FVC086) Size 3 64mm (FVC087) Size 4 70mm (FVC090) Size 5 76mm (FVC1593) Size 6 83mm (FVC098) Size 7 89mm (FVC100) Size 8 95mm (FVC103) Size 9 102mm (FVC105) Size 10 108mm (FVC106) Size 11 114mm (FVC1249) Size 12 121mm (FVC1250) Size 13 127mm (FVC1248) (sizes 11, 12 and 13 contain metal)
Description	Ring with support	Ring with support	Ring with support
Picture			
UOI	1	1	1
Route	Blue Diamond	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	More Significant Pelvic Organ Prolapse	More Significant Pelvic Organ Prolapse	More Significant Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	Pelvic Organ Prolapse	Mild 1 st degree Pelvic Organ Prolapse Cystocele	1 st degree Pelvic Organ Prolapse Mild Cystocele
Country of Manufacture	China	USA	USA


Ring with knob			
Supplier	RB Medical	Mediplus	Mediplus
NPC	Size 1 51mm (FVC1809) Size 2 57mm (FVC1624) Size 3 63mm (FVC1758) Size 4 70mm (FVC1608) Size 5 76mm (FVC1772) Size 6 83mm (FVC1734) Size 7 89mm (FVC1656)	Size 0 44mm (FVC314) Size 1 51mm (FVC315) Size 2 57mm (FVC316) Size 3 64mm (FVC317) Size 4 70mm (FVC318) Size 5 76mm (FVC319) Size 6 83mm (FVC320) Size 7 89mm (FVC321) Size 8 95mm (FVC322) Size 9 102mm (FVC323) Size 10 108mm (FVC1179) Size 11 114mm (FVC1180) Size 12 121mm (FVC1181) Size 13 127mm (FVC1178) (sizes 9, 10, 11, 12 and 13 contain metal)	Size 0 44mm (FVC1172) Size 1 51mm (FVC1169) Size 2 57mm (FVC1170) Size 3 64mm (FVC332) Size 4 70mm (FVC333) Size 5 76mm (FVC336) Size 6 83mm (FVC337) Size 7 89mm (FVC1904) Size 8 95mm (FVC1171) Size 9 102mm (FVC1177) Size 10 108mm (FVC1173) Size 11 114mm (FVC1174) Size 12 121mm (FVC1175) Size 13 127mm (FVC1176) (all sizes contain metal)
Description	Ring with knob	Ring with knob	Ring with knob
Picture			
UOI	1	1	1
Route	eDirect	Blue Diamond	Blue Diamond
Material	Flexible Silicone	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Any type of Pelvic Organ Prolapse + Stress Urinary Incontinence	Any type of Pelvic Organ Prolapse + Stress Urinary Incontinence	Any type of Pelvic Organ Prolapse + Stress Urinary Incontinence
Indications for Use (Manufacturer's)	Mild Pelvic Organ Prolapse w/ Mild Cystocele Urinary Incontinence	Stress Urinary Incontinence	Stress Urinary Incontinence
Country of Manufacture	USA	USA	USA

Ring with knob & support		
Supplier	RB Medical	Mediplus
NPC	Size 1 51mm (FVC1732) Size 2 57mm (FVC1793) Size 3 63mm (FVC1738) Size 4 70mm (FVC1798) Size 5 76mm (FVC1701) Size 6 83mm (FVC1711) Size 7 89mm (FVC1599)	Size 0 44mm (FVC1182) Size 1 51mm (FVC377) Size 2 57mm (FVC378) Size 3 64mm (FVC379) Size 4 70mm (FVC381) Size 5 76mm (FVC382) Size 6 83mm (FVC1245) Size 7 89mm (FVC384) Size 8 95mm (FVC1246) Size 9 102mm (FVC1247) Size 10 108mm (FVC1183) Size 11 114mm (FVC1242) Size 12 121mm (FVC1243) Size 13 127mm (FVC1244) (sizes 9, 10, 11, 12 and 13 contain metal)
Description	Ring with knob & support	Ring with knob & support
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Not included	Not included
Indications for Use (Manufacturer's)	Mild Pelvic Organ Prolapse w/ Mild Cystocele Urinary Incontinence	1 st -2 nd degree Pelvic Organ Prolapse Stress Urinary Incontinence Mild Cystocele
Country of Manufacture	USA	USA

Risser	
Supplier	Mediplus
NPC	Size 0 64x57mm (FVC1251) Size 1 70x62mm (FVC1252) Size 2 70x64mm (FVC1253) Size 3 60x67mm (FVC1254) Size 4 73x70mm (FVC1255) Size 5 78x79mm (FVC1256) Size 6 89x73mm (FVC1257) Size 7 89x84mm (FVC1258) Size 8 98x87mm (FVC1259) Size 9 105x86mm (FVC1260) (all sizes contain metal)
Description	Risser
Picture	
UOI	1
Route	Blue Diamond
Material	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	Stress Urinary Incontinence Mild Cystocele Cervical Incompetence Retrodisplacement
Country of Manufacture	USA

Shaatz		
Supplier	RB Medical	Mediplus
NPC	Size 0 38mm (FVC1730) Size 1 44mm (FVC1614) Size 2 51mm (FVC1659) Size 3 57mm (FVC1623) Size 4 64mm (FVC1633) Size 5 70mm (FVC1760) Size 6 76mm (FVC1755) Size 7 82mm (FVC1729) Size 8 88mm (FVC1672)	38mm (FVC1261) 44mm (FVC1262) 51mm (FVC1263) 57mm (FVC068) 64mm (FVC070) 70mm (FVC071) 76mm (FVC1264) 83mm (FVC075) 89mm (FVC1901) 95mm (FVC078)
Description	Shaatz with drain	Shaatz with drain
Picture		
UOI	1	1
Route	eDirect	Blue Diamond
Material	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Any type of Pelvic Organ Prolapse	Any type of Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	1 st -2 nd degree Pelvic Organ Prolapse Cystocele	1 st -2 nd degree Pelvic Organ Prolapse Vaginal Vault Mild Cystocele Rectocele
Country of Manufacture	USA	USA

Shelf			
Supplier	B Braun	Mediplus	Mediplus
NPC	Size 4 51mm (FVC1850) Size 5 54mm (FVC1823) Size 6 57mm (FVC1836) Size 7 60mm (FVC1822) Size 8 64mm (FVC1841) Size 9 67mm (FVC1842) Size 10 70mm (FVC1840) Size 11 73mm (FVC1824) Size 12 76mm (FVC1851) Size 13 79mm (FVC1849) Size 14 83mm (FVC1837) Size 15 86mm (FVC1815) Size 16 89mm (FVC1828) Size 17 92mm (FVC1827) Size 18 95mm (FVC1838) Size 19 98mm (FVC1826) Size 20 101mm (FVC1825)	Size 6 57mm (FVC1110) Size 8 64mm (FVC1111) Size 10 70mm (FVC1102) Size 12 76mm (FVC1104) Size 14 83mm (FVC1106) Size 16 89mm (FVC1109)	Size 6 57mm (FVC1189) Size 8 64mm (FVC1112) Size 10 70mm (FVC1103) Size 12 76mm (FVC1105) Size 14 83mm (FVC1107) Size 16 89mm (FVC1108)
Description	Shelf	Shelf with standard stem	Shelf with short stem
Picture			
UOI	1	1	1
Route	Blue Diamond	Blue Diamond	Blue Diamond
Material	Rigid Plastic	Flexible Silicone	Flexible Silicone
Indications for Use ¹	Advanced Pelvic Organ Prolapse	Advanced Pelvic Organ Prolapse	Advanced Pelvic Organ Prolapse
Indications for Use (Manufacturer's)	Pelvic Organ Prolapse Cystocele Rectocele	2 nd -3 rd degree Pelvic Organ Prolapse Procidentia	2 nd -3 rd degree Pelvic Organ Prolapse Procidentia
Country of Manufacture	UK	UK	UK

Smith	
Supplier	Mediplus
NPC	Size 0 64x57mm (FVC1251) Size 1 70x62mm (FVC1252) Size 2 70x64mm (FVC1253) Size 3 60x67mm (FVC1254) Size 4 73x70mm (FVC1255) Size 5 78x79mm (FVC1256) Size 6 89x73mm (FVC1257) Size 7 89x84mm (FVC1258) Size 8 98x87mm (FVC1259) Size 9 105x86mm (FVC1260) (all sizes contain metal)
Description	Smith (modified Hodge)
Picture	
UOI	1
Route	Blue Diamond
Material	Flexible Silicone
Indications for Use ¹	Refer to Manufacturer's Instructions
Indications for Use (Manufacturer's)	Stress Urinary Incontinence Mild Cystocele Incompetent Cervix Retrodisplacement
Country of Manufacture	USA