



May 1st 2024

Gloves Supply Disruption in Foley Catheter Trays

Dear Valued Customer,

This letter is to notify you that we are experiencing a supply disruption on the sterile gloves that come in our Foley catheter Trays. Due to the clinical need for these Trays, we will continue production without sterile gloves. The Trays will have a sticker on the outer packaging noting the lack of sterile gloves and notifying the nurse to obtain a sterile pair before beginning the catheterisation procedure.

Trays that do not contain the sterile gloves are expected to start entering the market around May 1st and the disruption is expected to last about 3 months. Impacted codes can be found in the table below. Please be assured we are doing everything we can to get sterile glove inventory from our suppliers so we can ensure you have complete catheter insertion trays for your patients.

We appreciate your patience as we work through this challenge.

Sincerely,

Rajesh Ramachandran
Market Access Manager – Urology and Critical Care

Table

Foley catheter	Collection system	Product code (Standard Length)
Biocath™ Hydrogel Coated Foley Catheter	Preconnected to 500mL Leg Bag	TR22645M (12–16)
	Preconnected to 2 Litre Bed Bag	TR22642L (12–16)
	Preconnected to Urinometer	TR2264UM (12–16)
Lubri-Sil™ Hydrogel Coated All Silicone Foley Catheter	Preconnected to 500mL Leg Bag	TR17585M (12–16)
	Preconnected to 2 Litre Bed Bag	TR17582L (12–16)
	Preconnected to Urinometer	TR1758UM (12–16)
Bardia™ PTFE Coated Latex Foley Catheter	Preconnected to 500mL Leg Bag	TR12455M (12–16)
	Preconnected to 2 Litre Bed Bag	TR12452L (12–16)
	Preconnected to Urinometer	TR1245UM (12–16)