

Information for Clinical Choice

Patient Temperature Management Blanket Warming Cabinets

Support Document

Introduction



Information for Clinical Choice (ICC) has been developed to assist Clinicians in the decision-making process when assessing the suitability of a product by providing a clear illustration and description of the features of a range of similar products supplied through NHS Supply Chain.

Blanket warming cabinets are used across hospitals including critical care areas, A&E resus rooms and theatre recovery, for warming blankets and any textiles for the maintenance of normothermia and the comfort of the patient.

The criteria provided, in the form of a Product Matrix and Support Document, is the result of a product review conducted by the clinical team in the Category Management Service Provider with support from clinical stakeholders from across the NHS as part of our assurance process.

The aim, alongside delivering savings back into NHS frontline services, is to ensure that clinical choice remains at the forefront of any product switching decision.

Patient Temperature Management: Blanket Warming Cabinets

| Blanket Warming Cabinets | | | | | | |
|--------------------------|---|---|---|---|---|---|
| Supplier | Central Medical Supplies Ltd | Central Medical Supplies Ltd | Central Medical Supplies Ltd | Central Medical Supplies Ltd | Central Medical Supplies Ltd | Central Medical Supplies Ltd |
| Brand | KanMed | Enthermics | Enthermics | Enthermics | KanMed | KanMed |
| MPC | GE2470 | EC250 | EC350 | EC750 | GE-2350S | GE-2380L |
| NPC | FSB85623 | FSB85628 | FSB85626 | FSB85625 | FSB85617 | FSB85605 |
| Description | Table Top Blanket cabinet | Blanket Warming cabinet 4-6 countertop | Blanket Warming Cabinet 10-12 blankets | Blanket Warming Cabinet 20-25 blankets | Warming cabinet Blanket with wheels 90cm | Blanket warming cabinet with wheels 180cm |
| Picture |  |  |  |  |  |  |
| Route | eDirect | eDirect | eDirect | eDirect | eDirect | eDirect |
| UOI | 1 | 1 | 1 | 1 | 1 | 1 |
| Temperature Range (Max) | 60°C | (71°C) | (71°C) | (71°C) | 30°C-80°C(80°C) | 30°C-80°C(80°C) |
| Country of Manufacture | Sweden | USA | USA | Sweden | Sweden | Sweden |

References:

If you have any questions, would like further information, or have feedback to share, please contact:

Raymond Pukacz - Clinical Engagement and Implementation Manager, Patient Temperature Management

raymond.pukacz@supplychain.nhs.uk

Keith Johnson - Category Manager, Patient Temperature Management

keith.johnson@supplychain.nhs.uk